

Nucleic Acid Purification

Genomic DNA purification kits - Tissue and cell

${\sf NucleoSpin}^{\circledR} {\sf Tissue\ procedure}$






Agarose gel electrophoresis of a 620bp amplified fragment from different

Genomic DNA was isolated from hair samples using NucleoSpin® Tissue 1: Marker 2: Gorilla gorilla graueri

- 3: Gorilla gorilla gorilla
- 4: Semnonitheous entellus
- 5: Trachypithecus poliocephalus (museum sample)



Analysis of genomic DNA from forensic samples (blood and epithelial cells of a rape victim, semen and blood of the suspect) Gel picture of amplified DNA from 3 STR loci, separated with the ABI 377 sequencer. Lane 15: Blood sample (victim) Lane 19: Fnithelial cells (victim) Lane 21: Epithelial/semen sample (victim/suspect) Lane 17: Blood sample (suspect) Sample preparation (differential lysis): After lysis of the epithelial cells and separation of semen cells the DNA of the suspect was purified with NucleoSnin® Tissue Data kindly provided by Dr S West, LKA Nordrhein Westfalen, Düsselforf,

DNA purification kits, genomic, NucleoSpin® Tissue



Total DNA extraction possible from: Tissue (e.g. mouse tails); Cells (e.g. bacteria); Clinical samples (stool, urine, biopsy samples); Forensic samples (dried blood spots, hair, buccal swabs, cigarette filters).

- Even very small amounts of DNA can be bound reversibly to the membrane (forensic analysis)
- Silica membrane technology
- Sample size: up to 25mg to tissue or 107 cells
- Yield: up to 35µg genomic DNA
- Elution volume: 60µL to 100µL
- Binding capacity: 60µg
- Preparation time: 20min/4 to 6 preps
- Format: mini spin columns
- Spin columns can be closed no cross-contamination
- No use of organic solvents
- Highly pure nucleic acids suitable for all common downstream applications

NucleoSpin® Tissue is designed for the rapid purification of highly pure genomic DNA from tissue samples, mouse tails, bacteria, yeast, forensic samples and clinical samples

Up to 35µg of high purity genomic DNA can be prepared (typical yields from tissue or cells: 15 to 25µg). The obtained DNA can be used directly for PCR*, Southern blotting or any kind of enzymatic reaction. With the NucleoSpin® Tissue method, lysis is achieved by incubation of the samples in a solution containing SDS and proteinase K at 56°C. Appropriate conditions for binding of DNA to the silica membrane of the NucleoSpin® Tissue columns are created by addition of large amounts of chaotropic ions (buffer B3) plus ethanol to the lysate. The binding process is reversible and specific for nucleic acids. Contaminations are removed by efficient washing. Pure genomic DNA is finally eluted under low ionic strength conditions in a slightly alkaline elution buffer and is ready to use for subsequent reactions.

Kit components: NucleoSpin® Tissue columns, 2mL collecting tubes, buffers, proteinase K.

For further information on this product please contact Customer Services, details can be found on the inside front cover.

Protocols are available for the isolation of Bacterial DNA, eg. Genomic DNA, eq. from Viral DNA, eq. CMV from stool Mycobacterium tuberculosis or Legionella Human or animal tissue pneumophila from sputum or brochoalveolar lavage CMV from urine EHEC bacteria from food (eg. fresh cows' milk) Mouse or rat tails Chlamydia trachomatis from cultures, biological fluids Bacteria or clinical specimens Yeast Borrelia burgdorferi from urine Dried blood spots (Guthrie cards) Hair roots Paraffin-embedded tissue Stool Insects Dental swabs

| Catalogue No | Alt. No | Description | Preps per kit |
|--------------|------------|------------------------|---------------|
| NZ74095210 | 740952.10 | NucleoSpin® Tissue kit | 10 |
| NZ74095250 | 740952.20 | NucleoSpin® Tissue kit | 50 |
| NZ740952250 | 740952.250 | NucleoSpin® Tissue kit | 250 |

Accessories

Buccal swabs

| Catalogue No | Alt. No | Description | Quantity |
|--------------|-----------|----------------------|----------|
| NZ74094025 | 740940.25 | Lysis buffer T1 25mL | 25mL |
| NZ740921 | 740921 | Wash buffer B5 20mL | 20mL |
| NZ740922 | 740922 | Wash buffer BW 100mL | 100mL |

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche