




## General Safety

### Disposal - glass



### Disposal bins, glass

02



Disposable, highly visible receptacle that helps reduce injuries by providing a safe laboratory work area. Made from sturdy corrugated cardboard, the receptacle has a 0.05mm thick polyethylene bag to contain broken glass fragments and spilled liquids. A pair of tabs keep the liner from slipping out of place. The snug-fitting lid has a flap through which waste material is easily inserted. When box is full, the internal safety cap is pulled into place and the entire unit is ready for disposal. The broken glass disposal box is printed with eye-catching words and graphics in four languages. Packaged flat to save space and shipping costs, it can be assembled in seconds.

Catalogue No	Alt. No	Description	Width, mm	Depth, mm	Height, mm	Pack qty
SAT-625-010R	F24653-0001	Floor type	300	300	680	6
SAT-625-020X	F24653-0002	Bench type	200	200	260	6



### Cartons, biohazard, incineration

02



Safe disposal and incineration of biohazard waste is simple and convenient with the Scienceware® biohazard incineration carton. The entire carton is designed to be burned along with its contents, thereby eliminating the need to transfer dangerous biohazardous material in plastic bags from work station to incinerator. This corrugated cardboard receptacle has a 0.04mm thick HMHD biohazard bag. A pair of tabs along the top edge of the carton prevent the bag from slipping out of place. The snug-fitting lid has a flap through which waste material is easily inserted. When the carton is filled, the internal safety cap is pulled into place and the entire unit is ready for incineration. Also features large, colourful biohazard graphics and safety instructions in four languages. Packed flat to save space and shipping costs, the carton may be assembled in seconds.

Catalogue No	Alt. No	Description	Width, mm	Depth, mm	Height, mm	Pack qty
SAT-627-010F	H13205-001	Floor type	300	300	600	6
SAT-627-020C	H13205-002	Bench type	200	200	250	6



### Disposal bin, aerosol/glass, Magpie™

04

Manufactured from strong corrugated board, for use with dry aerosol and glass. Specially made for landfill disposal. Letter box type opening. Flat packed.

Catalogue No	Capacity, L	Width, mm	Depth, mm	Height, mm	Pack qty
SAT-626-010L	10	275	275	450	10