

## Autoclave deodorisers, deodorant pads, Odo-Clave®



Impregnated with a heat released deodorant that neutralises autoclave odours.

Easy to handle and non messy.  
Pads are 32mm square.



Catalogue No	Legacy No	Fragrance	Pack qty
11976072	-	Clove	100
11956072	-	Cherry	100
11966072	-	Pine	100
12390469	-	Lemon	100

## Autoclave tape, indicator



For use with sterilisation packaging materials.

Meets BS 7893: 1997 specification for pressure sensitive adhesive closing and sealing tapes, also meets BS EN 867-2:1997 non-biological systems for use in sterilisers Part 2. Process indicators (Class A).



Catalogue No	Legacy No	Width, mm	Pack type
12340489	-	18	Roll/50m
12350489	-	24	Roll/50m

## Sterilisation indicators, Bowie and Dick Test, Comply™



For single use in porous load sterilisers for daily Bowie and Dick tests.

- Fast, efficient and disposable
- Simple and easy to read
- Stronger colours and higher definition
- Accurate, consistent, reliable and reproducible
- Easy record storage
- Saves time
- Paper can be re-used for note-taking



With unique colour chemistry, the 3M™ Comply™ Disposable Bowie and Dick type autoclave test pack conforms to BS EN 867-4 and carries the BSI KiteMark®.

Catalogue No	Legacy No	Description	Pack qty
12348199	-	Bowie and Dick test	20

## Sterilisation indicators, Bowie and Dick test



Single use Bowie and Dick type test pack. An accepted method of testing the steam penetration and air removal capability of a vacuum steriliser. Complies to BS EN ISO 11140-4:2007. BSI Kitemark® no. KM39542.



Catalogue No	Legacy No	Description	Pack qty
12380489	-	134°C to 138°C for 3.3 to 3.5min	20
12390489	-	121°C to 124°C for 8 to 8.3min	20

## Sterilisation indicators, TST™



Manufactured in accordance with ISO 11140 part 1, class 6. Class 6 specifies that the indicator shall have the capability to monitor the critical parameters of specified steriliser cycles. These cycles have been scientifically determined to give the best possible chance of achieving sterility.



Catalogue No	Legacy No	Description	Pack qty
12906996	-	121°C in 15min/134°C in 5.3min	100
12916996	-	134°C in 3.5min	100

Self adhesive

Catalogue No	Legacy No	Description	Pack qty
12926996	-	121°C in 15min/134°C in 5.3min (self adhesive)	100

## Paper, autoclave, Fisherbrand®



Green 60g/m² crepe paper for wrapping instruments prior to autoclaving.



Catalogue No	Legacy No	Dimensions, mm	Pack qty
11738513	-	500 x 500	500
11714586	-	1,000 x 1,000	250