

Bags, biohazard, polypropylene, autoclavable, roll packs, Fisherbrand®



Fisherbrand®

Standard temperature bags

- 40µm thickness
- With printed biohazard symbol
- Supplied on convenient rolls of 25 tear off bags per roll
- Suitable for temperatures of up to 121°C



Standard temperature (121°C), HDPE bags, rolls of 25 bags

Catalogue No	Legacy No	Width, mm	Depth, mm	Capacity, L	No. of rolls
11517762	-	250	400	3	40
11537762	-	300	500	9	30
11547762	-	400	780	21	16
11557762	-	500	800	33	20

High temperature (134°C) PP, bags in bulk packs

Catalogue No	Legacy No	Width, mm	Depth, mm	Capacity, L	Pack qty
11850482	-	400	780	21	500

Bags, biohazard, autoclavable, heavy duty, clear and printed



Scienceware®
BEL-ART PRODUCTS

Heavy duty bags that can be loaded with used pipettes, Petri dishes, culture plates and membrane filters.
The entire load is then sterilised at a maximum of 121°C and incinerated. All bags are 32µm thick.



Clear, unprinted can be used as conventional bin liners

Catalogue No	Legacy No	Width, mm	Height, mm	Pack qty
11392144	-	300	610	100
11302154	-	610	760	100
11312154	-	610	910	100

Clear, standard biohazard warning symbol and precautionary procedures in 4 languages (English, French, German and Spanish)

Catalogue No	Legacy No	Width, mm	Height, mm	Pack qty
11352154	-	300	610	100
11362154	-	610	760	100
11372154	-	610	910	100

Red, large black biohazard symbol and precautionary procedures in 4 languages (English, French, German and Spanish) with sterilisation indicator patch which turns black at 121°C

Catalogue No	Legacy No	Width, mm	Height, mm	Pack qty
11392154	-	360	480	200
11302164	-	480	580	200
11312164	-	640	890	200

Bags, biohazard, polyethylene, B-safe



Ideal for heavy duty biohazard disposal.

Supplied in handy, bench top dispenser boxes.

Catalogue No	Legacy No	Width, mm	Depth, mm	Pack qty
12300609	-	255	355	50
12310609	-	455	760	50
12320609	-	790	1,015	25