

Bottles PET



Bottles, PETG, diagnostic, sterile, Thermo Scientific Nalgene®



- Sterile
- Non-pyrogenic, non-cytotoxic

White lined closures are assembled onto each bottle. Bottles are tray packed, 20 per tray, 100 per case. Trays are packaged in double lined cartons for cleanroom use.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Closure size, mm	Pack qty
10672083	-	2035-0005	5	20	100
10704043	-	2035-0010	10	20	100
10682083	-	2035-0020	20	20	100



Bottles, PETG, transparent, square, narrow neck, sterile, Thermo Scientific Nalgene®



Radiation sterilised and non-pyrogenic. Reduced permeability to CO₂/O₂. Tamper evident seal for assurance of sterility. White HDPE screw closure. Packed in shrinkwrapped trays.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	Width, mm	Closure size, mm	Pack qty
10037950	-	2019-0030	30	14	63	38	20	96
10444481	-	2019-0060	60	18	82	41	24	96
10484481	-	2019-0125	125	28	108	52	38	48
10028010	-	2019-0250	250	28	147	58	38	48
11734414	-	2019-0500	500	28	177	73	38	24
11744414	-	2019-1000	1,000	28	218	93	38	24
10434491	-	2019-2000	2,000	39	270	116	95	12



Bottles, PET, sterile, octagonal



New polyethylene terephthalate (PET) bottles are ideal for storage and sampling of media, biologicals and other aqueous solutions.

- The PET resin makes the bottle lightweight, break resistant and as transparent as glass
- Packaged in two convenient, shrinkwrapped trays within an outer bag to assure cleanliness
- Break-resistant bottles are an ergonomic and easy to handle alternative to glass
- Bottles are sterilised to 10⁻⁶ SAL and certified non-pyrogenic
- Leakproof, tamper evident 31.7mm HDPE screw cap
- Moulded graduations for accurate measurements
- USP class VI compliant bottles are validated non-cytotoxic and non-hemolytic

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Height, mm	Closure style	Closure size, mm	Pack qty
11751595	-	431731	125	25	115	Screw cap	31.7	48
11761595	-	431732	250	50	142	Screw cap	31.7	48
11721595	-	431733	500	50	180	Screw cap	31.7	24
11741595	-	431734	1,000	50	220	Screw cap	31.7	24