

Bottles polypropylene



Bottles, polypropylene, translucent, wide neck, Thermo Scientific Nalgene®



For chemicals, specimens and general use.

- Autoclavable with cap loosened off threads
- Excellent for chemicals, materials and specimens where autoclaving is required
- Excellent durability and sample protection
- PPCO (Polypropylene co-polymer) offers better clarity than LDPE or HDPE
- Wide neck, easy to fill with dry materials or liquids
- Resistant to most acids, bases and alcohols
- Withstands temperatures up to +121°C

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10404671	-	2105-0001	30	23	64	36	28	12
11769813	-	2105-0002	60	20	86	38	28	12
10404501	-	2105-0004	125	31	99	51	38	12
10414501	-	2105-0008	250	33	130	61	43	12
10748331	-	2105-0016	500	43	168	74	53	12
10058240	-	2105-0032	1,000	53	201	91	63	6



Bottles, polypropylene, translucent, wide neck



Azlon

Wide neck bottles ideal for liquids and solids.

Manufactured from high clarity polypropylene with good chemical resistance. Clear graduations showing mL and fl oz. Large 'write-on' space for safe easy labelling, S.P.I. recycling codes moulded to base. Autoclavable to 121°C

Catalogue No	Legacy No	Alt. No	Capacity, mL	Height, mm	O.D., mm	Neck I.D., mm	Pack qty
11323494	-	BWP0060PGRAD	60	81	40	22	10
11333494	-	BLP0125PGRAD	125	104	50	32	10
11343494	-	BWP0250PGRAD	250	140	60	34	10
11353494	-	BWP0500PGRAD	500	161	75	44	10
11363494	-	BWP1000GRAD	1,000	205	90	53	5



Bottles, polypropylene, translucent, Mason jars, wide neck, Thermo Scientific Nalgene®



Can accept fixtures or closures with Mason bottle threads. Autoclavable.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10224501	-	2115-0500	500	64	175	75	70	6
10234501	-	2115-1000	1,000	64	211	94	70	6
10663062	-	2115-2000	2,000	64	241	122	70	4
10284251	-	2115-3000	3,000	64	246	155	70	1