

Burettes, borosilicate glass, single bore, interchangeable PTFE stopcock, blue ceramic graduations, Class A, BS-ISO 385: 2005, Fisherbrand®

Fisherbrand®

Catalogue No	Legacy No	Capacity, mL	Graduations, mL
10450893	-	25	0.10
10738081	-	50	0.10



Burettes, borosilicate glass, single bore, interchangeable PTFE stopcock, white enamel graduations, Class AS, BS-ISO 385: 2005



PYREX®

Complies with ISO 385, manufactured from precision bore borosilicate glass tubing for high accuracy. Fitted with self lubricating PTFE key with 1.5mm bore. Each burette has a +30s waiting time before final reading can be taken. With white enamel graduations for increased visibility with dark solutions.

Uncertified

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL
11397954	-	3295/12AS	10	0.02
11307964	-	3295/14AS	25	0.10
11317964	-	3295/16AS	50	0.10
11327964	-	3295/18AS	100	0.20

With certificate of calibration

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL
12463471	-	3295/22ASWC	10	0.02
12473471	-	3295/24ASWC	25	0.10
12483471	-	3295/26ASWC	50	0.10
12493471	-	3295/28ASWC	100	0.20



Service you expect

Our product range is only part of the story. You can be assured that you will receive service second to none...

- Large dynamic warehouse means high product availability
- State-of-the-art inventory technology ensures efficient order processing
- Customer Service team dedicated to making your contact with Fisher Scientific easy and efficient
- Knowledgeable field sales specialists to advise you in product choice
- Promotions and offers ensure you receive the best value for money

Contact your local Customer Service team

