

Tubes, centrifuge, Thermo Scientific Nalgene®



The polymethylpentene tubes are autoclavable and can withstand temperatures of up to 170°C. High chemical resistance. Cannot be used with naked flame. Can also be used as conventional test tubes. Crystal clear.

Polypropylene copolymer - without caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10058730	-	3110-0120	Round bottom	12	-	16 x 101 (approx.)	50,000	10
10623791	-	3110-0150	Round bottom	15	-	16 x 116 (approx.)	50,000	10
10460033	-	3110-0500	Round bottom	50	-	29 x 103 (approx.)	50,000	10
10152541	-	3110-1000	Round bottom	100	-	32 x 166 (approx.)	50,000	10
10374061	-	3103-0015	Conical	15	Yes	17 x 119 (approx.)	6,000	10
10602682	-	3103-0050	Conical	50	Yes	29 x 134 (approx.)	6,000	10

Oak Ridge tubes, polypropylene copolymer - with reusable, linerless non-contaminating 4mm screw caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10470241	-	3148-0050	Conical	35	5 to 35 x 1	28.5 x 113 (approx.)	50,000	10

Adapter - must be used with Thermo Scientific Nalgene 35mL conical Oak Ridge tubes

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10102401	-	DS3147-0050	Conical, for 10470241 and 10094690	-	-	-	-	6

Polycarbonate, conical - without caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10098690	-	3105-0015	Conical	15	Yes	17 x 119 (approx.)	6,000	10
10070701	-	3105-0050	Conical	50	Yes	29 x 134 (approx.)	6,000	10

Polycarbonate, round - without caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10614141	-	3117-0120	Round bottom	12	-	16 x 101 (approx.)	50,000	10
10088750	-	3117-0150	Round bottom	15	-	16 x 116 (approx.)	50,000	10
10150571	-	3117-9500	Round bottom	50	-	29 x 105 (approx.)	50,000	10
10557571	-	3117-0500	Round bottom	50	-	29 x 103 (approx.)	50,000	10

Polycarbonate bottles - with polypropylene linerless caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10589291	-	3122-0250	Bottle	250	-	82 x 125 (approx.)	27,500	4
10303971	-	3122-0500	Bottle	500	-	89 x 168 (approx.)	13,700	4
10313971	-	3122-1000	Bottle	1,000	-	98 x 177 (approx.)	7,100	4

Oak Ridge, polycarbonate - with polypropylene linerless caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10048840	-	3118-0010	Round bottom	10	-	16 x 80 (approx.)	50,000	10
10078670	-	3118-0030	Round bottom	30	-	26 x 92 (approx.)	50,000	10
10538341	-	3118-0050	Round bottom	43	-	29 x 114 (approx.)	50,000	10

Oak Ridge tubes, polycarbonate - with polypropylene caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10094690	-	3146-0050	Conical	35	Yes	28.5 x 104 (approx.)	50,000	10

entry continued