

Centrifuge tubes


Continuation

Tubes, centrifuge, Thermo Scientific Nalgene®



Polypropylene copolymer bottles - with caps

Catalogue No	Legacy No	Alt. No	Type	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
10343871	-	3120-0500	Bottle	500	-	73 x 169 (approx.)	4,800	4
10028850	-	3120-0250	Bottle	250	-	82 x 125 (approx.)	up to 13,200	4
10408312	-	3120-9500	Bottle	500	-	90 x 169 (approx.)	4,800	4
10242991	-	3120-1000	Bottle	1,000	-	98 x 174 (approx.)	7,100	4



Tubes, centrifuge

Kartell



Autoclavable. Polymethylpentene can withstand temperatures of up to 170°C. High chemical resistance. Cannot be used with naked flame. Can also be used as conventional test tubes. Crystal clear.

Polymethylpentene - without caps

Catalogue No	Legacy No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
12306307	-	1387	2	10	0.1	16 x 110 (approx.)	-	10
11371175	-	1388	2	15	0.2	18 x 117 (approx.)	-	10

Polypropylene - without caps

Catalogue No	Legacy No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
12336297	-	302	3	7	-	12 x 100 (approx.)	-	10
11361235	-	303	3	16	-	18 x 100 (approx.)	-	10
11371235	-	305	3	31	-	24 x 93 (approx.)	-	10
11341175	-	306	3	55	-	30 x 100 (approx.)	-	10
11351175	-	307	3	75	-	35 x 100 (approx.)	-	10
11361175	-	308	3	115	-	40 x 120 (approx.)	-	10
12396297	-	2300	2	10	0.1	16 x 110 (approx.)	-	10
11381175	-	300	2	10	-	16 x 108 (approx.)	1,300	10

Tubes, centrifuge

PYREX®



Pyrex® glass - without caps

Catalogue No	Legacy No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
11353544	-	1660/02M	1	10	No	16 x 90 (approx.)	3,000	1
11363544	-	1660/06M	1	25	No	24 x 90 (approx.)	3,000	1
11373544	-	3460/02M	2	15	Yes	17 x 110 (approx.)	3,000	1

Tubes, centrifuge



Conical. Double threaded caps 15mL and 50mL only.

Soda glass - without caps

Catalogue No	Legacy No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
12823462	-	-	2	13	-	17 x 111 (approx.)	1,300	10
12853472	-	-	3	15	-	17 x 111 (approx.)	1,300	10
12863472	-	-	2	15	0 to 10 x 0.1	17 x 111 (approx.)	1,300	10