

# Centrifuge tubes



## Tubes, centrifuge, 50mL, self-standing



CORNING

- Black printed graduations with large white marking area
- Certified non-pyrogenic, sterile and RNase/DNase free
- Sterile by gamma irradiation
- Bulk packed

PP = Polypropylene

Catalogue No	Legacy No	Alt. No	Material	Cap style	RCF, x g [max.]	Inner pack qty	Pack qty
10509891	-	430897	PP	Plug seal	3,000	25	500
10038980	-	430921	PP	Flat	3,000	25	500



## Tubes, centrifuge, 250mL/500mL



CORNING

- Conical tubes for larger volume centrifugation
- Each case contains a recyclable foam rack to facilitate handling
- Max. RCF 6,000xg
- Sterile and certified non-pyrogenic
- Support cushion must be used with this product unless the centrifuge rotor has appropriately shaped V-bottom holders

PP = Polypropylene

Catalogue No	Legacy No	Alt. No	Description	Material	Cap style	RCF, x g [max.]	Inner pack qty	Pack qty
10758961	-	430776	250mL tube	PP	Plug	6,000	6	102
10018910	-	431123	500mL tube	PP	Plug	6,000	6	36

### Accessories

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10202531	-	430236	Support cushion for 250mL tube	6
10789721	-	431124	Support cushion for 500mL tube	6



## Tubes, centrifuge, conical, 15mL, Thermo Scientific Sterilin



Thermo  
SCIENTIFIC

Manufactured in optically clear polystyrene or chemically resistant polypropylene, sterilised by gamma irradiation and certified non-pyrogenic to 0.25EU/mL. Leak tested to BS 5213. Recommended RCF values: polypropylene tubes 6,000xg, polystyrene tubes 3,500xg.

Catalogue No	Legacy No	Alt. No	Description	Pack qty
12087472	-	15PSB	Bulk wrapped, polystyrene	450
11369223	-	15PSR	Racked, polystyrene	450