

Tubes, microcentrifuge, polypropylene, Safe-Lock®



eppendorf

Designed for safe handling of valuable and toxic samples as well as radioactive substances or DNA.

- The small hook of the patented hinged lid clips around the rim of the tube giving trouble free heating
- Lid easily opened with one hand
- When open, tubes can be autoclaved at 121°C for 20 minutes
- Roughened writing surfaces (on side and lid)
- Biohazardous substances can be removed, free of aerosols, through a thin membrane in the centre of the lid with a cannula
- Autoclavable micropestle for resuspending pellets
- Available in Eppendorf® Biopur® guaranteed sterile, pyrogen-free, RNase-free, DNase-free, DNA-free, ATP-free
- Graduation scale
- Can be centrifuged at up to 23,500xg

The 0.5mL tube is optimised for use in PCR* experiments. Assorted colour packs are a mixture of yellow, red, blue, green and colourless.



*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche

Safe-Lock®

Catalogue No	Legacy No	Alt. No	Colour	Capacity, mL	Pack qty
10292611	-	0030121023	Colourless	0.5	500
10622731	-	0030121708	Assorted	0.5	500
10421193	-	0030121112	Yellow	0.5	500
10394591	-	0030121120	Red	0.5	500
10490993	-	0030121139	Blue	0.5	500
10768371	-	0030121147	Green	0.5	500
10458312	-	0030121155	Amber (light protection)	0.5	500
10509691	-	0030120086	Colourless	1.5	1,000
10548841	-	0030120694	Assorted	1.5	1,000
10708571	-	0030120159	Yellow	1.5	1,000
10282851	-	0030120167	Red	1.5	1,000
10767411	-	0030120175	Blue	1.5	1,000
10170381	-	0030120183	Green	1.5	1,000
10160431	-	0030120191	Amber (light protection)	1.5	1,000
10038760	-	0030120094	Colourless	2.0	1,000
10702132	-	0030120686	Assorted	2.0	1,000
10547381	-	0030120205	Yellow	2.0	1,000
10058770	-	0030120213	Red	2.0	1,000
10663981	-	0030120221	Blue	2.0	1,000
10430423	-	0030120230	Green	2.0	1,000
10252851	-	0030120248	Amber (light protection)	2.0	1,000

PCR* clean Safe-Lock®

Catalogue No	Legacy No	Alt. No	Colour	Capacity, mL	Pack qty
10318661	-	0030123301	Clear	0.5	500
10398031	-	0030123328	Clear	1.5	1,000
10073321	-	0030123344	Clear	2.0	1,000

Biopur® Safe-Lock® individually sealed

Catalogue No	Legacy No	Alt. No	Colour	Capacity, mL	Pack qty
10393431	-	0030121570	Natural	0.5	50
10518721	-	0030121589	Natural	1.5	100
10241931	-	0030121597	Natural	2.0	100

Accessory

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10683001	-	0030120973	Micropestle for 1.5mL and 2mL Safe-Lock® tubes and 1.7mL Safe-Twist screw cap tubes	10