

Centrifuge tubes

Tubes, microcentrifuge, polypropylene, snap cap, HydroLogix



Thermo
SCIENTIFIC

- All HydroLogix tubes are certified RNase/DNase, RNA/DNA, ATP and pyrogen free
- Clear
- Graduation markings
- Frosted flat cap and side

Catalogue No	Legacy No	Alt. No	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
11525564	-	3446	0.6	30,000	No	1,000
11559914	-	3449	0.6	30,000	Yes	1,000
11535564	-	3448	1.5	26,000	No	500
11569914	-	3451	1.5	26,000	Yes	500
11535534	-	3434	2.0	25,000	No	500
11579914	-	3453	2.0	25,000	Yes	500

Tubes, microcentrifuge, polypropylene, locking lid, HydroLogix



Thermo
SCIENTIFIC

- All HydroLogix tubes are certified RNase/DNase, RNA/DNA, ATP and pyrogen free
- Clear
- Graduation markings
- Frosted flat cap and side

These tubes feature an advanced locking system for safe and convenient centrifuging, boiling, freezing, incubation, or storage of samples for your most critical procedures.

The cap snaps closed, then the lock easily slips into place with one simple hand operation. The larger, rounder cap design provides more thumb contact while the smooth surface reduces the possibility of ripped or punctured gloves.

3455IW, 3457IW and 3459IW are individually wrapped.

Catalogue No	Legacy No	Alt. No	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
10021782	-	3454	0.6	30,000	No	1,000
12388202	-	3455-11	0.6	30,000	Yes	1,000
11838332	-	3455IW	0.6	30,000	Yes	500
11599914	-	3456	1.5	24,000	No	500
12398202	-	3457	1.5	24,000	Yes	500
11910937	-	3457IW	1.5	24,000	Yes	500
10656833	-	3458	2.0	25,000	No	500
11509924	-	3459	2.0	25,000	Yes	500
11920937	-	3459IW	2.0	25,500	Yes	500

For Thermo Scientific HyClone® media, please see the Fisher Scientific Life Science Catalogue.