

## Tubes, microcentrifuge, polypropylene, screw-cap, conical, HydroLogix



Suitable for vapour phase liquid nitrogen storage, boiling and autoclaving.

- Made from a crystal clear resin, your samples are easily identifiable when using these tubes
- Available in amber for light sensitive applications
- Conical bottoms allow for more complete sample retrieval
- Graduation markings
- The superior O-ring design forms a positive seal against the rim of the screw cap to guard against leakage
- All HydroLogix tubes are certified RNase/ DNase, RNA/DNA, ATP and pyrogen free
- Clear with frosted flat cap and side (3462, 3466 and 3470 are supplied without cap)

Catalogue No	Legacy No	Alt. No	Colour	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
11519924	-	3460	Clear	0.5	25,000	No	500
11930937	-	3460A	Amber	0.5	25,000	No	500
12378202	-	3431-11	Clear	0.5	25,000	Yes	500
11878342	-	3431A	Amber	0.5	25,000	Yes	500
13180553	-	3462	Clear	0.5	25,000	No	500
11549924	-	3464	Clear	1.5	20,000	No	500
11828342	-	3464A	Amber	1.5	20,000	No	500
11529924	-	3461	Clear	1.5	20,000	Yes	500
13191794	-	3461A	Amber	1.5	20,000	Yes	500
11569924	-	3466	Clear	1.5	25,000	No	500
11599924	-	3468	Clear	2.0	20,000	No	500
11509934	-	3468A	Amber	2.0	20,000	No	500
11539924	-	3463	Clear	2.0	20,000	Yes	500
13100563	-	3470	Clear	2.0	20,000	No	500



## Tubes, microcentrifuge, polypropylene, screw-cap, free standing, HydroLogix



Suitable for vapour phase liquid nitrogen storage, boiling and autoclaving.

- Made from a crystal clear resin, your samples are easily identifiable when using these tubes
- Available in amber for light sensitive applications
- Conical bottoms allow for more complete sample retrieval
- The superior O-ring design forms a positive seal against the rim of the screw cap to guard against leakage
- All HydroLogix tubes are certified RNase/ DNase, RNA/DNA, ATP and pyrogen free

With screw-cap and O-ring unattached

Catalogue No	Legacy No	Alt. No	Colour	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
11509944	-	3472	Clear	0.5	30,000	No	500
11519944	-	3472A	Amber	0.5	30,000	No	500
11529944	-	3474	Clear	1.5	18,000	No	500
11539944	-	3474A	Amber	1.5	18,000	No	500
11589944	-	3488	Clear	2.0	18,000	No	500
12521407	-	3488A	Amber	2.0	18,000	No	500

With screw-cap and O-ring attached

Catalogue No	Legacy No	Alt. No	Colour	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
11559924	-	3465	Clear	0.5	30,000	Yes	500
11838342	-	3465A	Amber	0.5	30,000	Yes	500
11589924	-	3467	Clear	1.5	18,000	Yes	500
11960937	-	3467A	Amber	1.5	18,000	Yes	500
11519934	-	3469-11	Clear	2.0	18,000	Yes	500
11868342	-	3469A	Amber	2.0	18,000	Yes	500

Without cap

Catalogue No	Legacy No	Alt. No	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
13110563	-	3422-11	0.5	30,000	No	500
13120563	-	3478	1.5	18,000	No	500



For Whatman VectaSpin™ filter centrifuge tubes please see Fisher Scientific Life Science Catalogue or visit our website

