

Tube rack, microcentrifuge, polypropylene, PCR* transfer tray, Fisherbrand®



Fisherbrand®

Save time when working with individual 0.2mL tubes or strips of 8 or 12 tubes.

Transfer tray is a thin strip of plastic designed to rest perfectly over a 96 well tube rack. Just use handles to lift and carry, ideal when using with thermal cyclers. Autoclavable. Tray measures 114mm x 74mm x 11mm. Supplied in a pack of assorted colours.



Catalogue No	Legacy No	Description	Pack qty
11737275	-	PCR* transfer tray, assorted colours (blue, green, purple, yellow and orange)	5

Tube racks, microcentrifuge, polypropylene, PCR* tube, reversible, Fisherbrand®



Fisherbrand®

Reversible polypropylene rack with removable hinged lid that snaps in place on either side of the rack.

PCR* side of the rack has 168 wells that hold 0.2mL tubes. The wells accommodate 8 or 12 tube strips or individual tubes. The opposite side of the rack has 40 wells that hold 0.5mL tubes. There are 12 wells for 1.5mL tubes on both sides of the rack. Autoclavable. Rack measures 206mm x 131mm x 54mm. Supplied in a pack of assorted colours.



Catalogue No	Legacy No	Description	Pack qty
11787265	-	Reversible PCR* rack, assorted colours (blue, green, purple, yellow and orange)	5

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffmann-La Roche

Tube racks, microcentrifuge, polypropylene, 96 well, reversible, Fisherbrand®



Fisherbrand®

Polypropylene racks have 96 wells to accommodate 0.5mL microtubes on one side, or 96 wells for 1.5mL to 2mL microtubes on the other side.

Arranged in standard 8 x 12 array. Handy lid snaps into place in the tabs on each end of rack. Grid lines are moulded on the side of the rack making it easier to see the wells. Raised rings are moulded on the opposite side of the rack making it easier to replace the 1.5mL to 2mL tubes from these wells. Autoclavable. Racks measure 246mm x 121mm x 50mm. Available in natural or assorted colours.



Catalogue No	Legacy No	Description	Pack qty
11728174	-	96 well reversible rack, assorted colours (blue, green, pink, yellow and orange)	5
11738174	-	96 well reversible rack, natural	5
11958014	-	96 well reversible rack, blue	5
11968014	-	96 well reversible rack, green	5
11978014	-	96 well reversible rack, yellow	5

Tube racks, floating, Fisherbrand®

Fisherbrand®

- Brightly coloured, HDR foam racks are ideal for floating tubes in water baths, tubs or beakers
- Closed cell foam resists water absorption so racks can be washed and reused
- Blue round rack features unique 'X' slits to hold a variety of tube sizes
- Include detachable carrying handles to lift samples easily



Catalogue No	Legacy No	Description	Colour
11704736	-	Holds 4 x 50mL tubes	Green
11714736	-	Holds 8 x 15mL tubes	Blue
11724736	-	Holds 24 x 1.5mL to 2.0mL tubes	Yellow
11736565	-	Holds 18 x 0.2mL to 2.0mL tubes	Blue