

Chromatography

vials and closures



Vials, shell, 1mL, 2mL and 4mL, Fisherbrand®

Fisherbrand®

Clear and amber Fisherbrand® shell vials are suitable for use with Alcott, Gilson, Shimadzu and Waters® (Wisp 96 and 48 position carousel) autosamplers.

Clear glass with transparent polyethylene plug

Catalogue No	Legacy No	Capacity, mL	Dimensions, mm	Pack qty
11581374	-	1.0	7.8 x 35	100
10145424	-	1.0	8 x 40	100
11555914	-	2.0	12 x 32	100
11516074	-	4.0	15 x 45	100

Amber glass with transparent polyethylene plug

Catalogue No	Legacy No	Capacity, mL	Dimensions, mm	Pack qty
10506075	-	1.0	8 x 40	100
11565914	-	2.0	12 x 32	100
12951201	-	4.0	15 x 45	100

Clear polypropylene with transparent polyethylene plug

Catalogue No	Legacy No	Capacity, mL	Dimensions, mm	Pack qty
11561404	-	1.0	40 x 8	100
12970941	-	4.0	44.6 x 14.65	100
12980941	-	3.0	44.6 x 14.65	100



Vials, crimp top and micro, 8mm, Fisherbrand®

Fisherbrand®

Fisherbrand® clear and amber glass crimp top and microvials are suitable for use with the following autosamplers: Agilent, Beckman, Carlo Erba, CTC, Fisons, Perkin Elmer, Shimadzu, Thermo and VWR (Merck®)/Hitachi, etc.

Clear glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11551364	-	Round bottom	0.3	5.5 x 31.5	100
11521374	-	Conical bottom	0.3	5.5 x 31.5	100
11571364	-	Conical bottom	0.7	7 x 40	100
11581364	-	Flat bottom	0.7	7 x 40	100
11561364	-	Flat bottom	0.8	8.2 x 30	100
11531374	-	Flat bottom	1.2	8.2 x 30	100

Amber glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11551374	-	Conical bottom	0.7	7 x 30	100
11541364	-	Conical bottom	0.7	7 x 40	100
11591364	-	Flat bottom	0.7	7 x 40	100
11541374	-	Flat bottom	1.2	8.2 x 40	100



Crimp seals, aluminium, 8mm, Fisherbrand®

Fisherbrand®

Catalogue No	Legacy No	Description	Thickness, mm	Pack qty
11501344	-	With virgin PTFE seal	0.25	100
11561344	-	With red-orange natural rubber/transparent TEF seal	1.00	100
11511344	-	With red-orange natural rubber/transparent TEF seal (Agilent quality)	1.00	100
10369304	-	With red rubber, beige seal (Agilent quality)	1.00	100
11571334	-	With red PTFE/white silicone/red PTFE seal	1.00	100
11591334	-	With white silicone/red PTFE seal	1.30	100
11551344	-	With white silicone/red PTFE seal with slit	1.30	100