

Chromatography

vials and closures



Vials, screw and micro, 8mm (small opening), Fisherbrand®

Fisherbrand®

Fisherbrand® screw top vials and microvials for use with the following autosamplers: Beckman, CTC, Gilson, Knauer, Shimadzu, Spark, Varian, VWR (Merck®)/Hitachi, etc.

Clear glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11565874	-	Vial, small opening	2.0	12 x 32	100
11525884	-	Vial, small opening, graduated with marking spot	2.0	12 x 32	100
11515884	-	Microvial, small opening	2.0	12 x 32	100

Amber glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11575874	-	Small opening	1.5	11.6 x 32	100
11595874	-	Small opening, graduated with marking spot	2.0	12 x 32	100



Micro inserts, for 2mL screw vials (small opening), Fisherbrand®

Fisherbrand®

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11762428	-	Clear glass, 15mm top	0.1	5 x 31	1,000
11861653	-	Clear glass, 9mm top	0.1	5 x 31	1,000
11858951	-	Clear glass, with assembled plastic spring	0.1	5 x 29	1,000
11772428	-	Clear glass, flat bottom	0.2	5 x 31	1,000



Caps, screw, polypropylene, 8mm, Fisherbrand®

Fisherbrand®

Centre hole, 5.5mm diameter

Catalogue No	Legacy No	Description	Thickness, mm	Pack qty
11511404	-	With red rubber/PTFE seal	1.0	100
11571384	-	With PTFE red/silicone white/PTFE red seal	1.0	100
11501394	-	With red-orange natural rubber/transparent TEF seal	1.3	100
11591394	-	With red butyl/grey PTFE seal	1.3	100
11561384	-	With white silicone/red PTFE seal	1.3	100
11571394	-	With white silicone/red PTFE seal, with slit	1.3	100

Pink caps, centre hole

Catalogue No	Legacy No	Description	Thickness, mm	Pack qty
12961011	-	Pink with red rubber/PTFE beige seal	1.0	100
12901021	-	Pink with PTFE red/silicone white/PTFE red seal	1.0	100
12911021	-	Pink with silicone white/PTFE red seal	1.3	100
12971011	-	Pink with silicone white/PTFE red seal	1.3	100
12981011	-	Pink with natural rubber red-orange/TEF transparent seal	1.3	100
12991011	-	Pink with butyl red/PTFE grey seal	1.3	100

Closed top

Catalogue No	Legacy No	Description	Thickness, mm	Pack qty
11511394	-	With red-orange natural rubber/transparent TEF seal	1.3	100
11501404	-	With red butyl/grey PTFE seal	1.3	100
11551394	-	With white silicone/red PTFE seal	1.3	100