

Vials, short thread, and microvials, 9mm (wide opening), Fisherbrand®



Due to their technical geometry, the vials can be used on all common autosamplers such as Agilent, HTA, Shimadzu, Thermo, Varian, Waters®, etc.



Clear glass

| Catalogue No | Legacy No | Description | Capacity, mL | Dimensions, mm | Pack qty |
|--------------|-----------|---|--------------|----------------|----------|
| 11575894 | - | Short thread vial with integrated 0.2mL micro insert, graduated with white area for labelling | 0.2 | 12 x 32 | 100 |
| 11575914 | - | TopSert® TPX short thread vial, with integrated 0.2mL micro insert | 0.2 | 12 x 32 | 100 |
| 12951011 | - | Short thread vial with integrated 0.3mL micro insert | 0.3 | 11.6 x 32 | 100 |
| 11535914 | - | Microliter vial | 0.9 | 11.6 x 32 | 100 |
| 11515894 | - | Microlitre vial, wide opening | 1.1 | 12 x 32 | 100 |
| 11565884 | - | Vial, wide opening | 2.0 | 12 x 32 | 100 |
| 11575884 | - | Vial, wide opening, graduated with white area for labelling | 2.0 | 12 x 32 | 100 |



Amber glass

| Catalogue No | Legacy No | Description | Capacity, mL | Dimensions, mm | Pack qty |
|--------------|-----------|--|--------------|----------------|----------|
| 10145714 | - | Vial with integrated 0.2mL micro insert with label and filling lines | 0.2 | 11.6 x 32 | 100 |
| 11515924 | - | TopSert™ TPX vial with integrated 0.2mL micro insert 15mm top | 0.2 | 11.6 x 32 | 100 |
| 11585884 | - | Vial, wide opening, graduated with white area for labelling | 1.5 | 12 x 32 | 100 |

Transparent polypropylene

| Catalogue No | Legacy No | Description | Capacity, mL | Dimensions, mm | Pack qty |
|--------------|-----------|-------------------------------|--------------|----------------|----------|
| 11565964 | - | Micro vial | 0.3 | 11.6 x 32 | 100 |
| 11585964 | - | TPX micro vial, crystal clear | 0.3 | 11.6 x 32 | 100 |
| 11545974 | - | Micro vial | 0.7 | 11.6 x 32 | 100 |
| 11505974 | - | Vial with filling lines | 1.5 | 11.6 x 32 | 100 |

Amber polypropylene

| Catalogue No | Legacy No | Description | Capacity, mL | Dimensions, mm | Pack qty |
|--------------|-----------|-------------------------|--------------|----------------|----------|
| 11515974 | - | Micro vial | 0.3 | 11.6 x 32 | 100 |
| 11535974 | - | Vial with filling lines | 1.5 | 11.6 x 32 | 100 |

Service you expect

Our product range is only part of the story. You can be assured that you will receive service second to none...

- Large dynamic warehouse means high product availability
- State-of-the-art inventory technology ensures efficient order processing
- Customer Service team dedicated to making your contact with Fisher Scientific easy and efficient
- Knowledgeable field sales specialists to advise you in product choice
- Promotions and offers ensure you receive the best value for money

Contact your local Customer Service team

