

# Chromatography

## vials and closures



### Vials, screw, 2mL, 10mm (wide opening), Fisherbrand®

**Fisherbrand®**


Fisherbrand® screw top vials are ideal for use with the following autosamplers: Jasco, Perkin Elmer, Shimadzu, Varian, Waters®, etc.

#### Clear glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11511474	-	Wide opening	2	12 x 32	100
11521474	-	Wide opening, graduated with marking spot	2	12 x 32	100

#### Amber glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11531474	-	Wide opening, graduated with marking spot	2	12 x 32	100



### Caps, screw, polypropylene, 10mm, with 7mm centre hole, Fisherbrand®

**Fisherbrand®**

Catalogue No	Legacy No	Description	Thickness, mm	Pack qty
11571474	-	With red PTFE/white silicone/red PTFE seal	1.0	100
11551474	-	With red-orange natural rubber/transparent TEF seal	1.3	100
11561474	-	With white silicone/red PTFE seal	1.3	100
11581474	-	With white silicone/blue PTFE seal, with cross slit	1.5	100
11501484	-	Closed top with natural rubber red-orange/TEF transparent seal	1.3	100



### Vials, screw, 4mL, 13mm opening, Fisherbrand®

**Fisherbrand®**

Fisherbrand® 4mL screw vials are ideal for use with the following autosamplers: Dionex, Shimadzu, Spark, Varian, VWR Merck®/Hitachi, Waters® (Wisp 48 position carousel), etc.

#### Clear glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11536044	-	Clear, 1st hydrolytic class glass	4	15 x 45	100
11576044	-	Clear glass, 1st hydrolytic class, graduated, with white area for marking	4	15 x 45	100

#### Amber glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11556044	-	Amber, 1st hydrolytic class glass	4	15 x 45	100
11586044	-	Amber glass, 1st hydrolytic class, graduated, with white area for marking	4	15 x 45	100