

Chromatography

vials and closures



Caps, screw, polypropylene, 24mm, UltraBond™, Fisherbrand®

Fisherbrand®

UltraBond™ caps and septa form an inseparable unit so that the septum liner cannot fall out or be pushed into the vial when being penetrated with a blunt needle. The 100% firm fit is achieved by a patented process without the use of any glue or adhesive but by changing the molecular structure of the surface of the cap and septa.

Catalogue No	Legacy No	Description	Septa material	Thickness, mm	Pack qty
12921221	-	White screw cap, 15mm centre hole	Natural silicone/beige PTFE (EPA quality)	3.2	100
12931221	-	White screw cap, closed top	Natural silicone/beige PTFE (EPA quality)	3.2	100
11530595	-	White screw cap, centre hole	Butyl red/PTFE grey	2.5	100
11540595	-	White screw cap, closed top	Butyl red/PTFE grey	2.5	100
11520595	-	White screw cap, centre hole	Silicone white/PTFE beige	3.2	100
12931001	-	White screw cap, closed top	Silicone natural/PTFE beige	3.2	100
12921001	-	White screw cap, closed top	PTFE/EPDM black/PTFE	2.0	100



Vials, precision thread, 18mm opening, headspace, Fisherbrand®

Fisherbrand®


Precision thread headspace vials are compatible with CTC PAL, Varian, Gerstel, Atas, Shimadzu and Agilent autosamplers.

Clear glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11506114	-	Precision thread headspace vial, clear glass, round bottom	10	22.5 x 46	100
12941221	-	Precision thread headspace vial, clear glass, round bottom	20	22.5 x 75.5	100

Amber glass

Catalogue No	Legacy No	Description	Capacity, mL	Dimensions, mm	Pack qty
11526114	-	Precision thread headspace vial, amber glass, round bottom	10	22.5 x 46	100
12951221	-	Precision thread headspace vial, amber glass, round bottom	20	22.5 x 75.5	100



Caps, screw, magnetic, precision thread, for headspace vials, for 18mm, Fisherbrand®

Fisherbrand®

8mm centre hole, silver

Catalogue No	Legacy No	Septa material	Thickness, mm	Pack qty
12961221	-	White silicone/red PTFE	1.3	100
12971221	-	Transparent blue silicone/white PTFE	1.3	100
12981221	-	White silicone/blue PTFE	1.5	100
12991221	-	Red butyl/grey PTFE	1.6	100
12980951	-	Silicone white/aluminium foil silver, open top	1.3	100
10607445	-	Silicone white/PTFE red, pre-cut star	1.5	100

Closed top, silver

Catalogue No	Legacy No	Septa material	Thickness, mm	Pack qty
12901231	-	White silicone/red PTFE	1.3	100
12911231	-	Red butyl/grey PTFE	1.6	100