

Chromatography

HPLC columns and accessories



HPLC Columns, Thermo Scientific Hypersil™, GOLD



Ideal for new method development and LC/MS applications.

- Exceptional peak symmetry and resolution
- Outstanding pH stability
- Based on next generation ultra pure silica
- Other sizes are available

C₁₈

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
10610564	-	25003-154630	4.6	3.0	150
10501695	-	25005-154630	4.6	5.0	150
10582465	-	25002-103030	3.0	1.9	100
10630204	-	25002-152130	2.1	1.9	150

C₁₈ Guard columns

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
11528270	-	25003-014001	4.0 and 4.6	3.0	10	4
11548270	-	25005-014001	4.0 and 4.6	5.0	10	4

C₈

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
10553035	-	25203-052130	2.1	3.0	50
10078952	-	25203-152130	2.1	3.0	150
10503235	-	25205-102130	2.1	5.0	100
10619405	-	25205-253030	3.0	5.0	250
10149494	-	25205-154030	4.0	5.0	150
10690564	-	25205-054630	4.6	5.0	50
10774414	-	25203-104630	4.6	3.0	100

C₈ Guard cartridges

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
10075814	-	25203-013001	3.0	3.0	10	4
10487694	-	25203-014001	4.0 and 4.6	3.0	10	4
10259644	-	25205-013001	5.0	3.0	10	4

C₄

Catalogue No	Legacy No	Alt. No	I.D., mm	Particle size, µm	Length, mm
10702717	-	25502-101030	1.0	1.9	100
10056074	-	25502-022130	2.1	1.9	20
10479054	-	25502-202130	2.1	1.9	200
10391585	-	25502-152130	2.1	1.9	150
10159534	-	25503-052130	2.1	3.0	50
10599735	-	25503-152130	2.1	3.0	150
10609585	-	25505-052130	2.1	5.0	50
10342585	-	25502-033030	3.0	1.9	30
10544105	-	25502-054630	4.6	1.9	50
10312635	-	25505-254630	4.6	5.0	250

C₄ Drop-in guard cartridges

Catalogue No	Legacy No	Alt. No	I.D., mm	Particle size, µm	Length, mm	Pack qty
10447534	-	25503-012101	2.1	3.0	10	4
10610626	-	25503-014001	4.6	3.0	10	4

aQ Columns

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
10621344	-	25302-101030	1.0	1.9	100
10582635	-	25302-102130	2.1	1.9	100
10776153	-	25305-052130	2.1	5.0	50
10111723	-	25302-053030	3.0	1.9	50
10074313	-	25303-103030	3.0	3.0	100
10611296	-	25305-253030	3.0	5.0	250
10497414	-	25305-153030	3.0	5.0	150
10629955	-	25303-054030	4.0	3.0	50
10008902	-	25303-154630	4.6	3.0	150
10150593	-	25305-104630	4.6	5.0	100
10725793	-	25303-034630	4.6	3.0	30
10582065	-	25305-254630	4.6	5.0	250