

HPLC Columns, Thermo Scientific, BETASIL™



- High surface area with high bonded phase coverage
- Wide range of chemistries on 100Å silica

C₁₈ Columns

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
10251275	-	70103-102130	2.1	3	100
10588965	-	70103-152130	2.1	3	150
10189604	-	70105-052130	2.1	5	50
10508015	-	70105-152130	2.1	5	150
10017204	-	70103-104630	4.6	3	100
10154783	-	70103-154630	4.6	3	150
10022063	-	70105-154630	4.6	5	150
10704746	-	70105-254630	4.6	5	250

C₁₈ Guard cartridges

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
10427334	-	70103-012101	2.1	3	10	4
10035513	-	70105-012101	2.1	5	10	4
10016214	-	70103-014001	4.0 to 4.6	3	10	4
10448094	-	70105-014001	4.0 to 4.6	5	10	4

C₈ Columns

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
10579365	-	70203-052130	2.1	3	50
10260955	-	70203-054630	4.6	3	50
10119184	-	70203-154630	4.6	3	150
10724746	-	70205-154630	4.6	5	150
10302775	-	70205-254630	4.6	5	250

C₈ Guard cartridges

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
10025684	-	70203-012101	2.1	3	10	4
10035684	-	70205-012101	2.1	5	10	4
10448864	-	70205-014001	4.0 to 4.6	5	10	4

Silica 100 columns

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
10198914	-	70005-052130	2.1	5	50
10557245	-	70005-154630	4.6	5	150
10057254	-	70005-254630	4.6	5	250

Silica 100 guard cartridges

Catalogue No	Legacy No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
10187644	-	70005-012101	2.1	5	10	4
10279174	-	70005-013001	3.0	5	10	4
10239224	-	70005-014001	4.0 to 4.6	5	10	4

