

Wipes, cleanroom, Class 10, ISO Class 4, Super Polx® 1200



- Continuous filament, non-run knitted polyester
- Laundered and packaged in a Class 10 cleanroom
- Extremely low in particle and fibre generation
- Low extractables
- Manufactured and tested with strict adherence to quality assurance guidelines
- Vacuum sealed




Catalogue No	Legacy No	Dimensions, mm	Pack qty
12807153	-	230 x 230	150
12817153	-	300 x 300	75

Contamination control mat, multi-layer, disposable, white



Low profile. Pack of 8 tack mats x 30 layers.



Catalogue No	Legacy No	Dimensions, mm	Pack qty
11555315	-	457 x 1,168 (approx.)	8

Face masks, cleanroom, Class 100, ISO Class 5, CR gamma



Face mask with apertured polyethylene film, available with ties or ear loops.

Gamma sterilised and individually packaged.



Catalogue No	Legacy No	Alt. No	Description	Pack qty
11715505	-	62467	With ties	200
10441891	-	62470	With ear loops	200

For a full range of safety products, please visit our website.