

Cylinders

borosilicate glass Class B



Cylinders, borosilicate glass, spouted, Class B, blue graduations

MBL®

DIN 12680, BS 604, ISO 4788

Easy to read blue graduations. Hexagonal base.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
12595686	-	CRC201	5	0.1	2
11945406	-	CRC202	10	0.2	2
11955406	-	CRC204	25	0.5	2
11965406	-	CRC206	50	1.0	2
11975406	-	CRC209	100	1.0	2
11985406	-	CRC211	250	2.0	2
11921473	-	CRC213	500	5.0	2
11990895	-	CRC214	1,000	10.0	2
11328064	-	CRC216	2,000	20.0	1



Cylinders, borosilicate glass, spouted, Class B, white graduations

PYREX®

DIN 12680, BS 604, ISO 4788

Manufactured from Pyrex® borosilicate glass, with white graduations for improved visibility with dark solutions. Hexagonal base.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
11373554	-	3200/02MP	5	0.1	2
11383554	-	3200/04MP	10	0.2	2
11393554	-	3200/06MP	25	0.5	2
11303564	-	3200/08MP	50	1.0	2
11313564	-	3200/10MP	100	1.0	2
10365380	-	3200/12MP	250	2.0	2
10203640	-	3200/14MP	500	5.0	2
11323564	-	3200/16MP	1,000	10.0	2
11333564	-	3200/18MP	2,000	20.0	1



Cylinders, borosilicate glass, spouted, Class B, white graduations

PYREX®

DIN 12680, BS 604, ISO 4788

Pyrex® borosilicate glass cylinders with detachable polyethylene base.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
10671821	-	3205/04MP	10	0.2	2
11343564	-	3205/06MP	25	0.5	2
11353564	-	3205/08P	50	1.0	2
11363564	-	3205/10P	100	1.0	2
11373564	-	3205/12P	250	2.0	2