

Electrochemistry

pH meters, electrodes and accessories

pH buffers, capsules, Thermo Scientific Russell



Each box of one value contains five hermetically sealed vials, each containing 10 buffer capsules. Each capsule makes 100mL buffer solution.

Catalogue No	Legacy No	Alt. No	Description	Pack qty
12963207	-	BC54	pH4	50
12973207	-	BC57	pH7	50
12983207	-	BC59	pH9	50
12993207	-	BC510	pH10	50



pH buffers, all in one kits, Thermo Scientific Orion®



For a wide variety of applications, and ideal for field testing due to its smaller size.

- Includes a bottle of pH electrode storage solution and pH4, pH7 and pH10 buffers
- With a bottle of non-hazardous general pH electrode cleaning solution

Catalogue No	Legacy No	Alt. No	Description
10432062	-	910199	Standard all in one pH buffer kit, 475mL each of pH4.01, pH7.00 and pH10.01 buffers, pH electrode storage solution, and pH electrode storage bottle
10669508	-	810199	All in one pH buffer kit, 475mL each of pH4.01, pH7.00 and pH10.01 buffers, pH electrode storage solution, 30mL pH electrode cleaning solution D, and pH electrode storage bottle
12973247	-	916099	All in one pH buffer kit, 60mL each of pH4.01, pH7.00 and pH10.01 buffers, 60mL pH electrode storage solution, and 60mL pH electrode cleaning solution D



Indicator paper, books, pH, Fisherbrand®



With reference chart on the book cover showing colour changes in steps of pH2.

Books of 20 leaves (50mm x 10mm).

Catalogue No	Legacy No	Range, pH	Pack qty
11532413	-	2 to 12	10



Indicator paper, books, litmus, Fisherbrand®



Books of 20 leaves (50mm x 10mm).

Catalogue No	Legacy No	Range, pH	Indicator	Colour change	Pack qty
11542413	-	5 to 8	Red	Red to blue	10
11552413	-	5 to 8	Blue	Blue to red	10
11562413	-	-	Neutral	To red or blue	10

