

Conductivity meter, bench top, SevenCompact™, S230



METTLER TOLEDO

Catalogue No	Legacy No	Alt. No	Description
11334233	-	30019033	Meter only: meter, electrode holder, protective cover, quick guide, declaration of conformity and test certificate
11314233	-	30019034	Meter kit: meter, electrode holder, protective cover, quick guide, declaration of conformity, test certificate with InLab® 731-ISM and 2 calibration sachets for 1,413µS/cm and 12.88mS/cm
11324233	-	30019035	USP/EP Kit: meter, electrode holder, protective cover, quick guide, declaration of conformity, test certificate with InLab® 741-ISM and 2 calibration sachets for 1,413µS/cm

Accessories

Catalogue No	Legacy No	Alt. No	Description
11334223	-	30014092	InLab® 731-ISM conductivity sensor, 4 pole graphite, epoxy shaft, 10µS/cm to 1,000mS/cm, fixed 1.2m cable, mini-DIN
11344223	-	30014093	InLab® 731-ISM conductivity sensor, 4 pole graphite, epoxy shaft, 10µS/cm to 1,000mS/cm, fixed 2m cable, Mini DIN
11354223	-	30014094	InLab® 741-ISM conductivity sensor, 2 pole VA4, stainless steel shaft, 0.01µS/cm to 500µS/cm, fixed 1.2m cable, Mini DIN
12946387	-	30014099	InLab® Trace kit, includes InLab® Trace conductivity sensor, 0.0001µS/cm to 1,000µS/cm, titanium shaft, Mini LTW and flow cell
10719335	-	11124300	RS-P25 printer with RS232 cable
11958011	-	51302876	LabX® direct pH for data transfer to PC sensor, 0.0001µS/cm to 1,000µS/cm, titanium shaft, mini-LTW and flow cell.
11384193	-	51191926	USB A-B cable (1.8m) for LabX® direct software
11394193	-	21901297	Barcode reader
11967951	-	21901309	USB cable for barcode reader
12820523	-	51302049	Conductivity standard, 1,413µS/cm, box of 30 x 20mL sachets
12764584	-	51302050	Conductivity standard, 12.88mS/cm, box of 30 x 20mL sachets



Have you thought about.....

Laboratory coat, unisex, Fisherbrand®



Manufactured to high quality standards, the Fisherbrand® unisex laboratory coat is exclusively available from Thermo Fisher Scientific. Manufactured using a quality pre-shrunk polycotton fabric, it is suitable for industrial laundering at temperatures up to 85°C.

Catalogue No	Legacy No	Size
11558102	-	S
11568102	-	M
11578102	-	L
11588102	-	XL
11598102	-	2XL
11508112	-	3XL



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.