

Filtration centrifugal tubes, Amicon® Ultra, 0.5mL



- High sample recovery (greater than 90% of dilute starting solution)
- Concentrates 500µL down to 15µL
- Concentration factor 25X to 30X
- Fast processing time - 10min to 30min
- Reverse spin capability provides consistent recoveries

With their vertical membrane design, Amicon® Ultra 0.5mL filters deliver great performance and lower spin times.

Technical Specification - General

Volume range, µL	15 to 20 (typical, final concentration)
Volume, µL	500 (maximum, initial sample)
RCF, x g	14,000 (concentration spin) 1,000 (reverse spin) (recommended)
RCF, x g [max.]	15,000
Membrane surface area, cm²	1 (active)
Hold up volume, µL	<1

Catalogue No	Legacy No	Alt. No	MWCO, kDa	Pack qty
10505274	-	UFC503008	30	8
10012584	-	UFC503024	30	24
10256744	-	UFC503096	30	96
10328054	-	UFC5030BK	30	500
10378434	-	UFC505008	50	8
10194304	-	UFC505024	50	24
10695445	-	UFC505096	50	96
10257324	-	UFC5050BK	50	500
10527164	-	UFC510008	100	8
10114694	-	UFC510024	100	24
11576515	-	UFC510096	100	96
10398244	-	UFC5100BK	100	500

Accessory

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10019193	-	UFC50VL96	Collection tubes	96



For the Fisherbrand® microcentrifuge, please see page 244



FINDING THE PRODUCTS YOU NEED

Have you thought about.....

Gloves, latex, disposable, Aloe Vera, powder free, Class I Medical, Fisherbrand®



Textured surface to provide ideal grip for dry and wet working conditions. Aloe Vera coating on the inner surface protects hands from drying, chapping and cracking, making this glove ideal for prolonged use.

Catalogue No	Legacy No	Size	Pack qty
11854873	-	S	100
11864873	-	M	100
11874873	-	L	100
11884873	-	XL	100



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.