

## Filtration media, quartz wool

Whatman®  
Part of GE Healthcare

Used for filtering acids and hot gases or as a plug for combustion tubes.

A mass of fibres, average diameter 6µm.

Catalogue No	Legacy No	Quantity, g
12373866	-	10.00



## Filtration media, floc, ashless grade

Whatman®  
Part of GE Healthcare

Assists coagulation by forming a more flocculent suspension.

Bulky mass of dry dispersed cellulose fibres. Ash content 0.015% by weight. Can also be used for forming a filter bed in a funnel.

Catalogue No	Legacy No	Alt. No	Quantity, g
10517121	-	1704-010	100.00



## Filtration media, paper clippings, ashless grade

Whatman®  
Part of GE Healthcare

Used as a filter bed for Buchner funnels for rapid filtration of coarse precipitates.

Offcuts of ashless filter papers which, when mixed with distilled water, form a suspension.

Catalogue No	Legacy No	Alt. No	Quantity, g
10624132	-	1703-050	500.00




## Have you thought about.....

### Cylinders, borosilicate glass, spouted, Class A, blue ceramic graduations, Fisherbrand®



Graduated in blue ceramic. Hexagonal base. Individual reorder code on each item.  
DIN 12680 BS 604 ISO 4788

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Pack qty
12952310	-	5	0.1	2
11517832	-	10	0.2	2
11527832	-	25	0.5	2
11537832	-	50	1.0	2
11547832	-	100	1.0	2
11557832	-	250	2.0	2



Look for the Fisherbrand® logo throughout the catalogue. It's your guide to great value products.