


Flasks, cell culture, polystyrene, Thermo Scientific Nunc®, Nunclon™ Δ treated



- Certified surface treatment for optimal cell attachment and growth
- Filter caps or vent/close caps available for good gas exchange
- Short, wide neck for easy access
- Individually leak tested
- Excellent optical quality
- Non-toxic, radiation sterilised
- 10205232 is barcoded with large code 128

Closure - vent/closed

Catalogue No	Legacy No	Alt. No	Culture area, cm ²	Neck type	Closure	Inner pack qty	Pack qty
10490033	-	163371	25	Angled	Vent/close	20	160
10068780	-	153732	80	Straight	Vent/close	5	50
10354451	-	156502	175	Straight	Vent/close	2	32
10568482	-	136196	25	Angled	Filter	20	160
10170961	-	178905	80	Straight	Filter	5	50
10127340	-	178883	175	Straight	Filter	2	32
10205232	-	178983	175	Straight	Filter	2	32

Carrier, disposable, for cell culture flasks, Fisherbrand®



- Lightweight and strong - holds four columns of up to five x 250mL cell culture flasks
- Rigid laminated cardboard
- Fits underneath a fume hood with generous openings for easy loading and unloading
- Printed area for sample identification
- Flat packed

Catalogue No	Legacy No	Description	Dimensions [w x d x h], mm	Pack qty
11938004	-	Disposable flask carrier, holds up to 20 x 75cm ² flasks	400 x 225 x 270	2

Choose online ordering

- Browse all our products and services with only a few clicks
- Take advantage of promotions throughout the year
- Orders are instantly and automatically processed
- Keep track of your budget using the order history functionality
- Access your contract prices
- Track your order online

UK and International: www.fisher.co.uk | **Ireland:** www.ie.fishersci.com

Sweden: www.fishersci.se | **Norway:** www.fishersci.no

Finland: www.fishersci.fi | **Denmark:** www.fishersci.dk

