



Funnel, for cell factory, polystyrene, Thermo Scientific Nunc®



- Pour cell suspension or media directly into cell factory
- Tubing not required
- Individually packed

Intended for filling cell factories with media or suspension in a LAF bench.

Catalogue No	Legacy No	Alt. No		Pack qty
10327192	-	140050		6



Disposable flasks, Erlenmeyer, PETG, sterile, Thermo Scientific Nalgene®



Optically clear PETG plastic. Ideal for shaker and suspension cell culture, media preparation and storage. Gamma sterilised, moulded-in graduations. Supplied with HDPE screw closures.

Plain bottom, white closure

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
10465622	-	4112-0125	125	25	24
10292032	-	4112-0250	250	50	12
10179661	-	4112-0500	500	50	12
10475622	-	4112-1000	1,000	50	6
10385012	-	4112-2000	2,000	500	4
10228732	-	4112-2800	2,800	250	4

Baffled bottom, white closure

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
10395012	-	4112-0125	125	25	24
10764033	-	4113-0250	250	50	12
10345112	-	4113-0500	500	50	12
10621513	-	4113-1000	1,000	50	6
10599592	-	4113-2000	2,000	500	4
10750664	-	4113-2800	2,800	250	4

Plain bottom, blue vented closure

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
10266432	-	4115-0125	125	25	24
10592033	-	4115-0250	250	50	12
10686673	-	4115-0500	500	50	12
10357302	-	4115-1000	1,000	50	6
10473582	-	4115-2000	2,000	500	4
10156192	-	4115-2800	2,800	250	4

Baffled bottom, blue vented closure

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
10529962	-	4116-0125	125	25	24
10403972	-	4116-0250	250	50	12
10082332	-	4116-0500	500	50	12
10206582	-	4116-1000	1,000	50	6
10434162	-	4116-2000	2,000	200	4
10228602	-	4116-2800	2,800	250	4

Vented closure, 0.2mm hydrophobic PTFE membrane

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10133532	-	4114-0038	Fits 250mL flasks	12