

Flasks

borosilicate glass, conical



Flasks, glass, conical, narrow neck, ISO 1773 and DIN 12353

PYREX®

Erlenmeyer, graduated.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
12579105	-	1130/03	5	12	52	10
12589105	-	1130/05M	10	12	60	10

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
11972568	-	1130/11M	25	18	70	10

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
11902619	-	1130/14D	250	34	145	10

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
11912619	-	1130/20D	500	34	180	10

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
12911921	-	1130/26D	1,000	42	220	10

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
12921921	-	1130/30D	2,000	50	280	10

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
12549115	-	1130/32D	3,000	52	310	2

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
10765332	-	1130/36D	5,000	52	365	1



Flasks, glass, conical, heavy duty, narrow neck

PYREX®

Erlenmeyer, graduated, specially designed to provide optimum mechanical strength. Reinforced rims to protect against accidental impact damage.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
12623072	-	4980-10	10	10	50	12
12580056	-	4980-25	25	13	65	12
12599995	-	4980-50	50	14	78	12
11941516	-	4980-125	125	23	114	12
11951516	-	4980-250	250	27	132	12

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
11931516	-	4980-500	500	31	176	6

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
11921516	-	4980-1L	1,000	37	216	6

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	Pack qty
12044777	-	4980-2L	2,000	43	268	1
10470160	-	4980-4L	4,000	43	360	1
10535200	-	4980-6L	6,000	43	410	1