

Flasks, borosilicate glass, conical, screw neck



PYREX®

Erlenmeyer, graduated, with rigid plastic screw cap. Ideal for the storage of culture specimens and media. The screw caps are fitted with an inert PTFE liner and rubber disc, and are resistant to repeated dry heat and steam sterilisation cycles, 121°C for 2hr.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Height, mm	Cap size, mm	Pack qty
11920905	-	1137/02M	50	85	22	10
11940905	-	1137/03D	100	105	30	10
11950905	-	1137/06M	250	130	42	10
11960905	-	1137/08D	500	160	42	10
10793032	-	1137/10M	1,000	210	42	1

Spare parts

Catalogue No	Legacy No	Alt. No	Description	Pack qty
12579835	-	4503/30	Caps for 100mL flasks	20
10181930	-	4503/42	Caps for 250mL, 500mL, and 1,000mL flasks	10



Flasks, borosilicate glass, Erlenmeyer, narrow neck, baffled




PYREX®

For use in the preparation of cultures.

Manufactured from Pyrex® borosilicate glass. With four baffled indents that provide optimal and consistent agitation.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Pack qty
11353454	-	1134/08	250	5
11363454	-	1134/10	500	5
11373454	-	1134/12	1,000	2
11383454	-	1134/14	2,000	1



Flasks, borosilicate glass, conical, wide neck, Fisherbrand®

Fisherbrand®

Erlenmeyer, graduated.

Catalogue No	Legacy No	Capacity, mL	Pack qty
12987620	-	50	10
12997620	-	100	10
12907630	-	250	10
12927630	-	500	10
11803452	-	1,000	1
11813452	-	2,000	1

