

Continuation

Incubators, orbital shaking, bench top, Thermo Scientific MaxQ™ HP420



Clamps

Catalogue No	Legacy No	Alt. No	Description	For use with
10474992	-	236010	Flask clamp, 25mL	MaxQ™ 8000 and HP series
10525433	-	236011	Flask clamp, 50mL	MaxQ™ 8000 and HP series
10729273	-	236012	Flask clamp, 125mL	MaxQ™ 8000 and HP series
10156272	-	236013	Flask clamp, 250mL/300mL	MaxQ™ 8000 and HP series
10535433	-	236014	Flask clamp, 500mL	MaxQ™ 8000 and HP series
10084032	-	236015	Flask clamp, 1L	MaxQ™ 8000 and HP series
10452902	-	236016	Flask clamp, 2L	MaxQ™ 8000 and HP series
10462902	-	236019	Fernbach flask clamp, 2,800mL	MaxQ™ 8000 and HP series
11826763	-	236025	Microplate clamp	MaxQ™ 8000 and HP series
11861331	-	236026	Utility tray, attaches to platform, 228.6mm x 457.2mm	MaxQ™ 8000 and HP series
11871331	-	236027	Utility tray, attaches to platform, 457.2mm x 457.2mm	MaxQ™ 8000 and HP series
10338922	-	236099	Flask clamp starter kit: 2 x 25mL, 2 x 50mL, 4 x 250/300mL, 4 x 500mL, 2 x 1L, 2 x 2L clamps and replacement screws	MaxQ™ 8000 and HP series
11736449	-	161033	Five tube inserts for 1L clamp, holds 3 x 50mL, 3 x 20mL, and 8 x 1.5L/2mL tubes	MaxQ™ 8000 and HP series
11706449	-	300349	Adhesive flask mat, 229mm x 229mm	MaxQ™ 8000 and HP series

Racks

Catalogue No	Legacy No	Alt. No	Description	For use with
12351653	-	600074	Test tube rack, 21mm to 25mm, 30mL tubes, 40 places	600088
10553503	-	600075	Test tube rack, 26mm to 30mm, 50mL tubes, 24 places	600088
11881331	-	600076	Test tube rack and holder for 10mm to 13mm tubes, 3mL to 7mL tubes, 72 places	MaxQ™ 8000 and HP series
10735544	-	600077	Test tube rack and holder for 16mm to 20mm tubes, 10mL to 20mL tubes, 40 places	MaxQ™ 8000 and HP series
11891331	-	600078	Test tube rack and holder for 21mm to 25mm tubes, 30mL tubes, 40 places	MaxQ™ 8000 and HP series
11769438	-	600079	Test tube rack and holder for 26mm to 30mm tubes, 50mL tubes, 24 places	MaxQ™ 8000 and HP series
10247562	-	600088	Universal adjustable angle test tube rack holder	MaxQ™ 8000 and HP series
10257562	-	600089	Microplate rack insert holds nine standard plates	MaxQ™ 8000 and HP series
10175752	-	600090	Microplate rack insert holds six deepwell plates	MaxQ™ 8000 and HP series
10503513	-	950040	Test tube rack, 10mm to 13mm, 3mL to 7mL tubes, 72 places	600088
10185612	-	950060	Test tube rack, 16mm to 20mm, 10mL to 20mL tubes, 40 places	600088

Accessories

Catalogue No	Legacy No	Alt. No	Description	For use with
10474432	-	452000	Dual stack shaker stand, includes kneewell side mount for orbital shaker stand	MaxQ™ HP 420 shakers
10545613	-	452001	Bench height single shaker stand	MaxQ™ HP 420 shakers
10664853	-	452002	Kneewell side mount, shaker slides out on rails for easy access	MaxQ™ HP 420 shakers
10719273	-	201127	Free standing 0 - 100°C, seven day chart recorder	MaxQ™ HP 420 shakers
13544940	-	197075	Replacement chart paper, 50 per box	MaxQ™ HP 420 shakers
12602795	-	194040	Quick change platform adapter	-
12612795	-	194140	Magnote wipe-off magnet and pen kit	Benchtop models only
13514940	-	194046	Clamp replacement screws	MaxQ™ 8000 and HP series
12652785	-	194045	Gassing manifold	MaxQ™ HP 420, 435/435/480/481