

Label printer, thermal transfer, BBP11



Ideal entry level thermal transfer printer for blood, tissue and slide labels in hospitals, healthcare and medical laboratories.

- Small and compact only 101.6mm (4inch) wide
- 300dpi print resolution
- Offers a wide range of printable materials
- Reliable and durable
- Many features built in as standard such as Ethernet (network), triple communication interface and dual walled construction

The Brady® BBP11 printer is a small and compact desktop printer which provides high quality thermal transfer printing at a low cost.

The printer offers high performance with durable and reliable output, with a printout speed of up to 127mm (5inch) per second. It is easy to use and offers quick ribbon and media loading. The standard internal Ethernet port enables easy integration into any application.

Catalogue No	Legacy No	Description
11583235	-	Benchtop thermal transfer label printer, BBP11

Accessories

Catalogue No	Legacy No	Description
10554335	-	76.2mm (3inch) external media carrier
10401125	-	R6400 ribbon, 65mm x 70m
10697865	-	R6400 ribbon, 110mm x 70m

FreezerBondz™ labels will adhere to already frozen samples.

Catalogue No	Legacy No	Description
12800159	-	FreezerBondz™, 25.4mm x 9.53mm



Labelling Software, LabelMark™



Compatible with Brady® printers BMP51, BMP53, BMP71, BBP11, IP300, laser printers and some other branded printers.

- Easy to use label creation software
- Create your own labels or import existing data
- Create 1D or 2D barcodes easily
- Suitable for use with Windows® 2000, XP, Vista and Windows® 7

Catalogue No	Legacy No	Alt. No	Description
10384460	-	LM 5	LabelMark™ software

