



Labels, microcentrifuge tube, Tough-Tags® and Tough-Spots®



Labels for microcentrifuge tubes and other containers in the research and clinical laboratory.

Unlike conventional tape materials, these labels withstand steam autoclaving, boiling water baths, freezer temperatures, organic solvents, caustic agents and other challenges without peeling.

Pre-cut, peel off labels perfectly sized to fit microcentrifuge tubes and other containers in the research and clinical laboratory. Chemically inert polyester with strong adhesion to plastic and other materials. Microcentrifuge tubes with Tough-Tags® wrap around labels slide easily in and out of rotors without binding because labels are thin and the adhesive does not migrate even after boiling and autoclaving. Available in laser sheets for labelling of large columns of samples direct from a laser printer.

Tough-Tags®, rolls

Catalogue No	Legacy No	Description	Colour	Dimensions, mm	Pack qty
12338807	-	Small	Clear	25 x 12.7	1,000
12348807	-	Small	White	25 x 12.7	1,000
12378807	-	Small	Red	25 x 12.7	1,000
12388807	-	Small	Yellow	25 x 12.7	1,000
12398807	-	Small	Green	25 x 12.7	1,000
12308817	-	Small	Blue	25 x 12.7	1,000
12358807	-	Large	Clear	32.5 x 12.7	1,000
12368807	-	Large	White	32.5 x 12.7	1,000
12318817	-	Large	Red	32.5 x 12.7	1,000
12328817	-	Large	Yellow	32.5 x 12.7	1,000
12338817	-	Large	Green	32.5 x 12.7	1,000
12348817	-	Large	Blue	32.5 x 12.7	1,000

Tough-Spots®, rolls

Catalogue No	Legacy No	Description	Colour	Dimensions, mm	Pack qty
12358817	-	Small	White	9.5	1,000
12368817	-	Small	Red	9.5	1,000
12378817	-	Small	Yellow	9.5	1,000
12388817	-	Small	Green	9.5	1,000
12398817	-	Small	Blue	9.5	1,000
12308827	-	Large	White	12.7	1,000
12318827	-	Large	Red	12.7	1,000
12328827	-	Large	Yellow	12.7	1,000
12338827	-	Large	Green	12.7	1,000
12348827	-	Large	Blue	12.7	1,000

Cryogenic labels for laser printers, sheets

Catalogue No	Legacy No	Description	Colour	Dimensions, mm	Pack qty
12672580	-	Tough Spots™ large	White	9.5	3,840
12662580	-	Tough-Tags™ small	White	24 x 12.7	2,975
12682580	-	Tough-Tags™ large	White	32.5 x 12.7	2,125



Label dispenser, Write-On™



With convenient writing platform. The epoxy coated platform's finger holes and built in cutting edge make dispensing and cutting the tape a simple one handed action.

Includes a plastic core for 76mm I.D. tape rolls up to 25.4mm wide. Tapes not included.

Catalogue No	Legacy No	Description
11397154	-	Label dispenser, for Write-On™ label tape