

Microscope slides, glass ring cells



For mounting thick specimens such as foraminifera, insect wings, etc. on microscope slides. May also be used for observing Brownian movement.

Catalogue No	Legacy No	Diameter, mm	Pack qty
12322158	-	15	10



See Lens Cleaning
Tissue page 1022




Chamber slides, TC-treated, Thermo Scientific Nunc™, Lab-Tek™, Chamber Slide™



Developed for microscopy of in-situ tissue culture.

- Cells grow on a microscope slide
- Growth chamber can be removed when culturing is complete
- Useful for viral and mycoplasma testing, chromosome studies, toxicity tests and immunocytology
- Certified surface treatment ensures optimal conditions for cell attachment and growth
- Permanox™ and glass do not interfere with fluorescent labels

TPX = Permanox™



Catalogue No	Legacy No	Alt. No	No. of wells	Slide material	Volume, working, mL	Culture area, cm ²	Inner pack qty	Pack qty
13043043	-	177372	1	Glass	2.0 to 4.5	9.4/well	8	96
11377764	-	177410	1	TPX	2.5 to 4.5	9.4/well	8	96
11347764	-	177380	2	Glass	1.2 to 2.0	4.2/well	8	96
11636141	-	177380SP	2	Glass	1.2 to 2.0	4.2/well	8	16
11387764	-	177429	2	TPX	1.2 to 2.0	4.2/well	8	96
11357764	-	177399	4	Glass	0.5 to 0.9	1.8/well	8	96
12832794	-	177399SP	4	Glass	0.5 to 0.9	1.8/well	8	16
11397764	-	177437	4	TPX	0.5 to 0.9	1.8/well	8	96
11367764	-	177402	8	Glass	0.2 to 0.4	0.8/well	8	96
13053043	-	177402SP	8	Glass	0.2 to 0.4	0.8/well	8	16
11307774	-	177445	8	TPX	0.2 to 0.4	0.8/well	8	96
11327774	-	178599	16	Glass	0.1 to 0.2	0.4/well	8	96
11317774	-	177453	Flaskette	Glass	2.5 to 5.0	10/well	8	96
10439672	-	170920	SlideFlask	PS	2.5 to 5.0	9.0/well	5	50
12862794	-	177445SP	8	TPX	0.02 to 0.4	0.8/well	8	16

Accessories

Catalogue No	Legacy No	Alt. No	Description	Inner pack qty	Pack qty
11387734	-	171862	Cover slip for 10439672	100	1,000
11357734	-	171080	Cover slip for 11327774	-	96

All products are for research use only, unless otherwise indicated.