


Microscope slide mailers, Fisherbrand®

Fisherbrand®

Range of slide mailers to hold standard slides of 75mm x 25mm, 2 or 5 place in natural coloured polypropylene.

Catalogue No	Legacy No	Description	Dimensions, mm	Pack qty
11755088	-	2 place polypropylene, colour natural	93 x 68 x 5	25
11719885	-	5 place polypropylene, end opening	43 x 24 x 88	25



Microscope slide dispenser, Fisherbrand®

Fisherbrand®

Easy to load and dispense 75mm x 25mm standard slides. Turning either side wheel dispenses slides one at a time and dispenser holds up to 91 slides.

Catalogue No	Legacy No	Dimensions, mm
11761556	-	113 x 123 x 163



Microscope slide mailers

Kartell

Particularly useful for cytological and geological specimens.

Robust protection of 76mm x 26mm specimen slides during transportation. Each slide has its own compartment with support bars at 55mm centre to protect subject matter from being scuffed. 'Top' is clearly marked for ease of opening. Constructed of polypropylene.

Catalogue No	Legacy No	Alt. No	No. of places	Width, mm	Depth, mm	Height, mm
11328623	-	923	1	50	100	5
11341185	-	925	2	71	84	5
11338623	-	924	3	84	100	5



Storage tray rack, microscope slide folders, Fisherbrand®



Fisherbrand®

- Unique rack allows the easy storage and retrieval of five 20 place slide folders (11703217)
- Stackable with folder identification clasps for reference and quick retrieval
- Ideal for room temperature or refrigerator storage of samples

Technical Specification - General

Temperature range, °C	-10 to 90
Dimensions [w x d x h], mm	226 x 340 x 125

Catalogue No	Legacy No	Description
11720594	-	Storage tray rack for microscope slide folders



All products are for research use only, unless otherwise indicated.