

Pipettes, bulb, soda lime glass, BS 1583 Class A, Fisherbrand®

Fisherbrand®

Colour coded to BS 3996.

| Catalogue No | Legacy No | Capacity, mL | Tolerance, mL | Colour | Pack qty |
|--------------|-----------|--------------|---------------|--------|----------|
| 11536943 | - | 1 | ±0.008 | Blue | 2 |
| 11546943 | - | 2 | ±0.01 | Orange | 2 |
| 11556943 | - | 5 | ±0.015 | White | 2 |
| 11566943 | - | 10 | ±0.02 | Red | 2 |
| 11576943 | - | 25 | ±0.03 | Blue | 2 |
| 11586943 | - | 50 | ±0.05 | Red | 2 |



Pipettes, bulb, soda lime glass, one mark, ISO 648, Class AS



MBL®

- Available with or without works certificate
- Complies with ISO 648 and DIN 12691
- Class AS calibrated for delivery with +5s waiting time
- Printed with batch code for downloadable batch certificate – works certified pipettes have each pipette coded with individual serial number
- Works certified pipettes supplied with calibration certificate which allows traceability to national standards
- Markings in permanent amber stain
- Manufactured from soda lime glass with tooled jets for extra strength and to ensure correct outflow times

Without certificate of calibration

| Catalogue No | Legacy No | Alt. No | Capacity, mL | Tolerance, mL | Colour | Pack qty |
|--------------|-----------|----------|--------------|---------------|--------|----------|
| 11369643 | - | PRA720AS | 1 | ±0.008 | Blue | 2 |
| 11379643 | - | PRA722AS | 2 | ±0.010 | Orange | 2 |
| 11389643 | - | PRA724AS | 5 | ±0.015 | White | 2 |
| 11399643 | - | PRA726AS | 10 | ±0.020 | Red | 2 |
| 11309653 | - | PRA728AS | 15 | ±0.030 | Green | 2 |
| 11319653 | - | PRA730AS | 20 | ±0.030 | Yellow | 2 |
| 11329653 | - | PRA732AS | 25 | ±0.030 | Blue | 2 |
| 11339653 | - | PRA734AS | 50 | ±0.050 | Red | 2 |
| 11349653 | - | PRA736AS | 100 | ±0.080 | Yellow | 2 |

With certificate of calibration

| Catalogue No | Legacy No | Alt. No | Capacity, mL | Tolerance, mL | Colour | Pack qty |
|--------------|-----------|------------|--------------|---------------|--------|----------|
| 12464681 | - | PRA740ASWC | 1 | ±0.008 | Blue | 2 |
| 12484681 | - | PRA742ASWC | 2 | ±0.010 | Orange | 2 |
| 12404691 | - | PRA744ASWC | 5 | ±0.015 | White | 2 |
| 12424691 | - | PRA746ASWC | 10 | ±0.020 | Red | 2 |
| 12444691 | - | PRA748ASWC | 15 | ±0.030 | Green | 2 |
| 12464691 | - | PRA750ASWC | 20 | ±0.030 | Yellow | 2 |
| 12484691 | - | PRA752ASWC | 25 | ±0.030 | Blue | 2 |
| 12404701 | - | PRA754ASWC | 50 | ±0.050 | Red | 2 |
| 12424701 | - | PRA756ASWC | 100 | ±0.080 | Yellow | 2 |

