


Pipettes, straight, serological, glass, sterile



- Graduated glass pipettes

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Length, mm	Jet	Inner pack qty	Pack qty
12314448	-	D804/PS	1	0.01	230	Short	10	500
12324448	-	D805/PS	1	0.10	230	Short	10	500
11778245	-	D802/PS	1	0.01	270	Drawn	10	500
12344448	-	D806/PS	2	0.02	260	Short	10	500
11798245	-	D803/PS	2	0.02	270	Drawn	10	500
12364448	-	D807/PS	5	0.05	310	Drawn	5	500
10377864	-	D808/PS	10	0.10	310	Drawn	5	250

Pipettes, straight, serological, polystyrene, Fisherbrand®



Suitable for tissue culture work.

A black magnifier stripe runs lengthwise along the pipette making the meniscus obvious and reducing reading errors. The pipettes have excellent optical clarity with high contrast black markings making it easier to read volumes. The 5mL, 10mL, 25mL and 50mL sizes have both ascending and descending graduations. Tip design on the 25mL and 50mL sizes increases the effect of surface tension to minimise dripping. The individually wrapped pipettes come in a clear, all plastic wrap and the multipacks come in a recyclable LDPE bag.



Individually wrapped

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
11819660	-	1	0.01	Yellow	100	1,000
11567742	-	2	0.01	Green	100	500
11577742	-	5	0.1	Blue	50	200
11839660	-	10	0.1	Orange	50	200
11517752	-	25	0.2	Red	50	200
11537752	-	50	0.5	Purple	25	100

Multipack

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Colour	Inner pack qty	Pack qty
11889650	-	1	0.01	Yellow	25	1,000
11869650	-	2	0.01	Green	25	500
11899650	-	5	0.1	Blue	50	500
11809660	-	10	0.1	Orange	25	500
12559546	-	25	0.2	Red	20	200
11859660	-	50	0.5	Purple	20	100