

# Pipettes straight



## Pipettes, straight, serological, polystyrene, paper/plastic wrap, Stripette®



CORNING

Ideal for sterile tissue culture applications, especially while wearing latex gloves as the wrap reduces static cling.

- Individually wrapped
- A strip creates a thermometer like effect at the meniscus, making the volume easy to read
- Colour coding on wrapping allows easy size identification
- Sterile, non-pyrogenic and DNase and RNase free
- Exclusive anti drip tip assures accurate delivery
- Bidirectional graduations provide choice of ascending and descending scales

| Catalogue No | Legacy No | Alt. No | Volume, mL | Colour        | Inner pack qty | Pack qty |
|--------------|-----------|---------|------------|---------------|----------------|----------|
| 10074450     | -         | 4485    | 1          | Yellow stripe | 50             | 1,000    |
| 10288370     | -         | 4486    | 2          | Green stripe  | 50             | 1,000    |
| 10420201     | -         | 4487    | 5          | Blue stripe   | 50             | 200      |
| 10084450     | -         | 4488    | 10         | Orange stripe | 50             | 200      |
| 10606151     | -         | 4489    | 25         | Red stripe    | 25             | 200      |
| 10041591     | -         | 4490    | 50         | Purple stripe | 25             | 100      |
| 10255221     | -         | 4491    | 100        | Aqua stripe   | 10             | 100      |



## Pipettes, straight, serological, polystyrene, plastic wrap, bulk pack, Stripette®



CORNING

- Sterile, non-pyrogenic and DNase/RNase free
- Exclusive anti-drip tip assures accurate delivery
- Bidirectional graduations provide choice of ascending and descending scales

| Catalogue No | Legacy No | Alt. No | Volume, mL | Colour        | Inner pack qty | Pack qty |
|--------------|-----------|---------|------------|---------------|----------------|----------|
| 10420672     | -         | 4010    | 1          | Yellow stripe | 50             | 1,000    |
| 10574712     | -         | 4020    | 2          | Green stripe  | 50             | 1,000    |
| 10758782     | -         | 4050    | 5          | Blue stripe   | 50             | 500      |
| 10275471     | -         | 4100    | 10         | Orange stripe | 50             | 500      |
| 10633053     | -         | 4250    | 25         | Red stripe    | 50             | 200      |
| 10534961     | -         | 4500    | 50         | Purple stripe | 50             | 100      |



## Pipettes, straight, serological, polystyrene, plastic wrap, Stripette®



CORNING

Ideal for sterile applications such as tissue culture.

- Feature a unique vertical strip which creates a thermometer like effect at the meniscus, making volume easy to read
- 'Banana peel' plastic wrap provides easy, aseptic opening and pipette removal
- Individual, clear, plastic wrap allows full view of pipette for easy size identification and reduces storage space requirements
- 10194331 has a wide tip for handling viscous fluids
- Sterile, non-pyrogenic and DNase and RNase free
- Exclusive anti drip tip assures accurate delivery
- Bidirectional graduations provide choice of ascending and descending scales

| Catalogue No | Legacy No | Alt. No | Volume, mL | Colour        | Inner pack qty | Pack qty |
|--------------|-----------|---------|------------|---------------|----------------|----------|
| 10238800     | -         | 4011    | 1          | Yellow stripe | 100            | 1,000    |
| 10265471     | -         | 4012    | 1          | Yellow stripe | 100            | 200      |
| 10340501     | -         | 4021    | 2          | Green stripe  | 100            | 1,000    |
| 10127400     | -         | 4051    | 5          | Blue stripe   | 50             | 200      |
| 10677341     | -         | 4101    | 10         | Orange stripe | 50             | 200      |
| 10206391     | -         | 4492    | 10         | Orange stripe | 50             | 200      |
| 10732742     | -         | 4251    | 25         | Red stripe    | 50             | 200      |
| 10636391     | -         | 4501    | 50         | Purple stripe | 25             | 100      |
| 10194331     | -         | 4484    | 100        | Aqua          | 10             | 100      |

