

## Pipette fillers, Thermo Scientific Matrix S1



- Separate aspirate and dispense controls
- Lithium ion battery technology
- Backlit LCD display
- 1mL pipette stabiliser
- Customisable ID area

Includes table stand, wall mount, 1mL pipette stabiliser, universal charging transformer.

### Technical Specification - General

Aspirate speed, s	8
Dispensing speed, s	8 (+ gravity dispense)
Battery	Lithium ion
Run time, hr	20 (continuous use between charges approx.)
Recharge time, min	60 (approx)
Weight, g	220
Nosepiece Pipette Holder	Autoclavable silicone rubber
Pipette Compatibility	All major brands of glass or plastic serological pipettes, 1mL to 100mL
Filter assembly, µm	0.45
Fill time, s	5 (for 50mL pipette)



Catalogue No	Legacy No	Alt. No	Colour
10402822	-	9501	White
10729534	-	9511	Clear
10626293	-	9521	Blue
10367302	-	9531	Red
10072332	-	9541	Green

### Accessories

Catalogue No	Legacy No	Alt. No	Description
12964527	-	9057	Matrix CM, hydrophobic filters, 0.45µm, pack of 24
10030524	-	9064	Matrix CM2 nosepiece holder, silicone gripper, filter
10200124	-	9065-11	Matrix CM pipette gripper, silicone, pack of 4
11681950	-	9066	Matrix CM2 wall mount
10335524	-	9067	Matrix CM table stand
11866244	-	9068	Matrix CM2 power supply, universal, 100-240V 50/60Hz
10108553	-	9070	Matrix CM2 pipette support, 1mL

## Pipette fillers, PipetBoy Acu



**INTEGRA**

Cable free, pipette controller.

Extremely lightweight, compact and easy to use. Works on a rechargeable 9V NiMH battery which is easily inserted at the top of the handle. One charge is sufficient to fill and purge at least 500 x 50mL pipettes. Accepts standard pipettes ranging from 0.1mL to 100mL. Comes in a range of five transparent colours and a traditional light grey to aid identification. With battery, transformer for 230V/50Hz, sterile filter and wall mount.

Catalogue No	Legacy No	Alt. No	Colour
10013451	-	155000	Light grey
10195551	-	155015	Green
10165601	-	155016	Red
10798252	-	155017	Blue
10768082	-	155018	Translucent
10060262	-	155019	Purple

### Accessories

Catalogue No	Legacy No	Alt. No	Description
10339051	-	153005	Rechargeable battery 9V NiMH
10574982	-	151005	Handle mount
10541941	-	151020	Pipettor, rubber, 0.1mL to 50mL
10531941	-	153015	Filter, sterile, for PipetBoy