

Pipettor tips, filter, universal fit, Fisherbrand® SureOne™




Fisherbrand® SureOne™ filter pipettor tips are a comprehensive line of universal fit pipettor tips, available in a range of volumes, from 5µL to 5mL, in sterile racked and sterile reload packaging. Compatible with most popular brands of pipettor, SureOne™ achieves optimal fit with Fisherbrand® Elite™ pipettors. The SureOne™ filter tip reload system offers the same benefits as the SureOne™ standard pipettor tip reload system and is the only stackable filter tip reload system, providing an environmental solution without compromising sterility.

Reload system


- Offset wells allow filter tips to nest and stack
- Rigid inserts offer a robust non-pliable design providing secure tip seal
- Lot numbered inserts for traceability
- Easy to read filter tip volume
- Sterile packaging with outer sleeve protects tips from exterior contamination
- Environmentally friendly filter tip packaging option reduces plastic waste and saves bench top space
- Blister sleeve guide rails align tip inserts for smooth transfer into the rack

Non-sterile products - certified to be free of RNase/DNase and DNA.

Sterile products - e-beam sterile products are certified to be free of RNase/DNase, pyrogen, bioburden, PCR* inhibitors and endotoxins.



Racked sterile



5µL nano micropoint tip, graduated at 2.5µL

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12932561	-	0.1 to 5	Filtered sterile	Clear	96	960
12942561	-	0.1 to 5	Filtered sterile reload	Clear	96	960



10µL micropoint tip, graduated at 2µL

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12922541	-	0.1 to 10	Filtered sterile	Clear	96	960
12952561	-	0.1 to 10	Filtered sterile reload	Clear	96	960



10µL extended length micropoint tip, graduated at 2.5µL

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
11512813	-	0.1 to 10	Filtered sterile	Clear	96	960
12962561	-	0.1 to 10	Filtered sterile reload	Clear	96	960



20µL extended length micropoint tip, graduated at 2µL, 5µL, 10µL and 20µL (use 12992531 10µL to 20µL empty rack box for reloads)

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12972561	-	0.1 to 20	Filtered sterile	Clear	96	960
12982561	-	0.1 to 20	Filtered sterile reload	Clear	96	960



20µL universal beveled tip

(use 12992541 20µL to 300µL empty rack box for reloads)

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
11552813	-	2 to 20	Filtered sterile	Clear	96	960
12992561	-	2 to 20	Filtered sterile reload	Clear	96	960



100µL universal beveled tip

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12992541	-	10 to 100	Filtered sterile	Clear	96	960
12902571	-	10 to 100	Filtered sterile reload	Clear	96	960



entry continued

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffmann-La Roche