

Pipettors, variable volume, Thermo Scientific Finnpiptette® F2



- Large ergovisio display with white numbers on black background for increased visibility
- Double action pipetting button, rotating upper part prevents accidental volume alteration
- Light, smooth pipetting action enables better results over longer periods
- Fully autoclavable, entire pipettor can be autoclaved
- 'Super Blow Out' piston with volumes of 50µL and below ensures delivery of micro drops

Technical Specification - Specific

Colour code	10073453 Pink	10686614 Pink	11827351 Pink	10596264 Yellow	10050464 Turquoise	10413865 Yellow	12326542 Turquoise
Volume range, µL	0.2 to 2 micro	0.5 to 5 micro	1 to 10 micro	1 to 10	2 to 20 micro	2 to 20	5 to 50 micro
Increments, µL	0.002	0.01	0.02	0.02	0.02	0.02	0.02
Accuracy, %	±12.0 to ±2.5	±6.0 to ±1.5	±2.5 to ±1.0	±3.5 to ±1.0	±3.0 to ±1.0	±3.0 to ±1.0	±3.0 to ±0.6
Precision, %	10.0 to 2.0	5.0 to 1.0	2.0 to 0.5	3.0 to 0.8	2.5 to 0.4	2.5 to 0.4	2.5 to 0.3
Colour code	12336542 Yellow	10554474 Yellow	10155433 Yellow	10541834 Blue	10762984 Green	10329263 Red	
Volume range, µL	5 to 50	10 to 100	20 to 200	100 to 1,000	500 to 5,000	1,000 to 10,000	
Increments, µL	0.02	0.2	0.2	1	10	20	
Accuracy, %	±3.0 to ±0.6	±3.0 to ±0.8	±1.8 to ±0.6	±1.0 to ±0.6	±2.0 to ±0.5	±2.0 to ±0.5	
Precision, %	2.5 to 0.3	1.0 to 0.2	0.7 to 0.2	0.6 to 0.2	0.8 to 0.2	0.8 to 0.2	

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Finntip®	Colour code
10073453	-	4642010	0.2 to 2 micro	0.002	Flex 10, 10, 50	Pink
10686614	-	4642020	0.5 to 5 micro	0.01	Flex 10, 10, 50	Pink
11827351	-	4642030	1 to 10 micro	0.02	Flex 10, 10, 50	Pink
10596264	-	4642040	1 to 10	0.02	Flex 200, 250 Univ	Yellow
10050464	-	4642050	2 to 20 micro	0.02	50 micro	Turquoise
10413865	-	4642060	2 to 20	0.02	Flex 200, 250 Univ	Yellow
12326542	-	4642120	5 to 50 micro	0.02	50 micro	Turquoise
12336542	-	4642130	5 to 50	0.02	Flex 200, 300, Univ 250, 200 Ext	Yellow
10554474	-	4642070	10 to 100	0.2	Flex 200, 250 Univ, 300, 200 Ext	Yellow
10155433	-	4642080	20 to 200	0.2	Flex 200, 250 Univ, 300, 200 Ext	Yellow
10541834	-	4642090	100 to 1,000	1	Flex 1000, 1000, 1000 Ext	Blue
10762984	-	4642100	500 to 5,000	10	5mL	Green
10329263	-	4642110	1,000 to 10,000	20	10mL, Flex 10mL Ext	Red

Accessory

Catalogue No	Legacy No	Alt. No	Description
10693195	-	9420400	F-stand, linear stand, white, 6 positions

