

Pipettors, variable volume, multichannel, Thermo Scientific Finnpiquette® F2



- Large ergovisio display with white numbers on black background for increased visibility
- Double action pipetting button, rotating upper part prevents accidental volume alteration
- Light, smooth pipetting action enables better results over longer periods
- Fully autoclavable, entire pipettor can be autoclaved
- 'Super blow out' piston with volumes of 50µL and below ensures delivery of micro drops

Technical Specification - Specific

	10228063	10115483	10125483	10339263	10500865
No. of channels	8	8	8	8	12
Volume range, µL	1 to 10	5 to 50	10 to 100	30 to 300	1 to 10
Accuracy, %	±12.0 to ±2.4	±5.0 to ±1.5	±5.0 to ±1.3	±5.0 to ±1.0	±12.0 to ±2.4
Precision, %	8.0 to 1.6	2.0 to 0.7	2.0 to 0.5	2.0 to 0.3	8.0 to 1.6

	10305854	10315854	10794545	10401444	10598554
No. of channels	12	12	12	16	16
Volume range, µL	5 to 50	10 to 100	30 to 300	1 to 10	5 to 50
Accuracy, %	±5.0 to ±1.5	±5.0 to ±1.3	±5.0 to ±1.0	±12.0 to ±2.4	±5.0 to ±1.5
Precision, %	2.0 to 0.7	2.0 to 0.5	2.0 to 0.3	8.0 to 1.6	2.0 to 0.7



8 channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Finntip®	Colour code
10228063	-	4662000	1 to 10	0.02	Flex 10, 10 micro, 50 micro	Pink
10115483	-	4662010	5 to 50	0.1	Flex 200, 250 Univ, 200 Ext	Yellow
10125483	-	4662020	10 to 100	0.2	Flex 200, 250 Univ, 200 Ext	Yellow
10339263	-	4662030	30 to 300	1.0	Flex 300, 300	Orange

12 channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Finntip®	Colour code
10500865	-	4662040	1 to 10	0.02	Flex 10, 10 micro, 50 micro	Pink
10305854	-	4662050	5 to 50	0.1	Flex 200, 250 Univ, 200 Ext	Yellow
10315854	-	4662060	10 to 100	0.2	Flex 200, 250 Univ, 200 Ext	Yellow
10794545	-	4662070	30 to 300	1.0	Flex 300, 300	Orange

16 channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Finntip®	Colour code
10401444	-	4662080	1 to 10	0.02	20 micro, 50 micro	Purple
10598554	-	4662090	5 to 50	0.1	50 micro	Turquoise

Accessory

Catalogue No	Legacy No	Alt. No	Description
10693195	-	9420400	F-stand, linear stand, white, 6 positions