

Pipettors, variable volume, single and multichannel, mLINE®



Multichannel

	13137853	13157504	13107903
No. of channels	8	8	8
Volume range, µL	0.5 to 10	5 to 100	30 to 300
Volume, µL	10 5 1	100 50 10	300 150 30
Inaccuracy, %	1.50 2.50 4.00	0.70 1.00 3.00	0.60 1.00 2.00
Imprecision, %	1.00 2.50 4.00	0.25 0.70 1.50	0.25 0.50 1.00

	13187933	13157554	13184424
No. of channels	12	12	12
Volume range, µL	0.5 to 10	5 to 100	30 to 300
Volume, µL	10 5 1	100 50 10	300 150 30
Inaccuracy, %	1.50 2.50 4.00	0.70 1.00 3.00	0.60 1.00 2.00
Imprecision, %	1.00 2.50 4.00	0.25 0.70 1.50	0.25 0.50 1.00

Single channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Tip, I.D.
13127284	-	725010	0.1 to 3	0.002	A, L, M
13181564	-	725020	0.5 to 10	0.01	A, L, M
13117254	-	725030	2 to 20	0.02	C, N
13104754	-	725050	10 to 100	0.10	C, D, O, V
13197963	-	725060	20 to 200	0.20	C, D, P, Q
13137973	-	725070	100 to 1,000	1.00	E, G, S
13151674	-	725080	500 to 5,000	10.0	J
13164414	-	725090	1,000 to 10,000	20.0	U

8 channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Tip, I.D.
13137853	-	725120	0.5 to 10	0.01	A, L, M
13157504	-	725130	5 to 100	0.10	C, D, O, V
13107903	-	725140	30 to 300	0.20	D, Q

12 channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Tip, I.D.
13187933	-	725220	0.5 to 10	0.01	A, L, M
13157554	-	725230	5 to 100	0.10	C, D, O, V
13184424	-	725240	30 to 300	0.20	D, Q

Accessories

Catalogue No	Legacy No	Alt. No	Description
13127923	-	725600	Carousel stand for 6 mechanical pipettors
13194424	-	725620	Linear stand for all Sartorius Biohit pipettors

