

Pipettors, electronic, variable volume, eLINE®



12 channel

Catalogue No	Legacy No	Alt. No	Volume range, µL	Increments, µL	Tip, I.D.
13171874	-	730423	0.2 to 10	0.05	A, L, M
13128093	-	730443	5 to 120	0.50	C, D, V
13197284	-	730463	10 to 300	1.00	D, Q
13117903	-	730493	50 to 1,200	5.00	H, T, Z

Accessories

Catalogue No	Legacy No	Description
13157544	-	Charging carousel for four pipettors
13104254	-	Charging stand for one pipettor
13101754	-	Charging adapter, 230V, 50Hz
12479780	-	Battery pack, NiMH



Pipettor tips, standard, Optifit



Sartorius Biohit tips manufactured from natural coloured, virgin, non-wettable polypropylene. Supplied in single trays, stacked trays or bulk packs. Pre-sterilised tips are shrink wrapped and gamma irradiated. CE IVD.

Catalogue No	Legacy No	Alt. No	Volume range, µL	Tip, I.D.	Pack type	Inner pack qty	Pack qty
13154234	-	790010	0.1 to 10	A	Single trays	96	960
13167594	-	790011	0.1 to 10	A	Single trays sterile	96	960
12449750	-	790012	0.1 to 10	A	Refill tower	96	960
13138023	-	790013	0.1 to 10	A	Refill pack sterile	96	1,920
13167204	-	790014	0.1 to 10	A	Bulk in box	-	1,000
13127274	-	790200	0.5 to 200	C	Single trays	96	960
13161684	-	790201	0.5 to 200	C	Single trays sterile	96	960
13101804	-	790202	0.5 to 200	C	Refill tower	96	960
13154784	-	790203	0.5 to 200	C	Refill pack sterile	96	1,440
13154304	-	790204	0.5 to 200	C	Bulk in box	-	1,000
13184714	-	790350	5 to 350	D	Single trays	96	960
13117973	-	790351	5 to 350	D	Single trays sterile	96	960
13184554	-	790352	5 to 350	D	Refill tower	96	960
13174554	-	790354	5 to 350	D	Bulk in box	-	1,000
13127973	-	791000	10 to 1,000	E	Single trays	96	960
13117923	-	791001	10 to 1,000	E	Single trays sterile	96	960
13174744	-	791002	10 to 1,000	E	Refill tower	96	960
13134424	-	791005	10 to 1,000	E	Bulk in box	-	400
13151884	-	771004	10 to 1,000	E	Bulk in box	-	1,000
13164274	-	791020	10 to 1,000	G	Single trays	96	960
13164684	-	791021	10 to 1,000	G	Single trays sterile	96	960
12929637	-	791024	10 to 1,000	G	Bulk in box	-	1,000
13147524	-	791200	50 to 1,200	H	Single trays	96	960
13104414	-	791201	50 to 1,200	H	Single trays sterile	96	960
12489760	-	791202	50 to 1,200	H	Refill tower	96	960
13197614	-	791203	50 to 1,200	H	Refill pack sterile	96	960
13157264	-	791204	50 to 1,200	H	Bulk in box	-	1,000
13164504	-	791210	50 to 1,200	Z	Single trays	96	960
13184524	-	791211	50 to 1,200	Z	Single trays sterile	96	960
13194584	-	791212	50 to 1,200	Z	Refill pack	96	960
13107624	-	791213	50 to 1,200	Z	Refill pack sterile	96	960
12419770	-	780304	100 to 5,000	J	Single trays	50	500
12999647	-	780305	100 to 5,000	J	Single trays sterile	50	500
13197304	-	780300	100 to 5,000	J	Bulk in box	-	100
13164454	-	780308	100 to 5,000	J	Bulk in box	-	1,000
13121874	-	780310	0.5mL to 10mL	K	Bulk in box	-	100

Accessories

Catalogue No	Legacy No	Alt. No	Description	Pack qty
13137284	-	790910	Tip box for refill system for 10µL, 300µL and 350µL tips (Tip I.D. A, C and D)	10
13107254	-	720920	Tip box for refill system for 1,000µL and 1,200µL tips (Tip I.D. E and H)	10

