

Pipettors, fixed volume, single or two stroke, Ovation®



Fixed volume models are the ideal ergonomic choice for settings where pipetting volumes never deviate from an exact protocol.

- Ergonomic design
- Filtered nozzle to prevent liquid contamination
- Requires less force to discard tips than traditional pipettors
- Low contamination risk, the pipettor stands securely on the counter
- Both single stroke and two stroke (overflow) models are available in a number of volume increments ranging from 1µL to 1,000µL

Single stroke

Catalogue No	Legacy No	Alt. No	Colour	Volume, µL	Accuracy within ±%	Precision, %
11912225	-	1157-0005	Orange	5	1.5	0.60
11922225	-	1157-0010	Orange	10	1.5	0.50
11932225	-	1157-0020	Yellow	20	1.0	0.30
11942225	-	1157-0025	Green	25	1.4	0.60
11952225	-	1157-0030	Green	50	0.8	0.24
11962225	-	1157-0100	Green	100	0.8	0.15
11972225	-	1157-0200	Blue	200	0.8	0.15
11982225	-	1157-0250	Purple	250	1.2	0.30
11992225	-	1157-0300	Purple	300	1.0	0.25
11902235	-	1157-0500	Purple	500	0.8	0.20
11922235	-	1157-1000	Purple	1,000	0.8	0.15

Two stroke

Catalogue No	Legacy No	Alt. No	Colour	Volume, µL	Accuracy within ±%	Precision, %
11912245	-	1257-0001	Red	1	2.7	1.80
11922245	-	1257-0002	Red	2	1.5	0.90
11932245	-	1257-0005	Orange	5	1.5	0.60
11956885	-	1257-0010	Orange	10	1.5	0.50
11966885	-	1257-0020	Yellow	20	1.0	0.30
11976885	-	1257-0025	Green	25	1.4	0.60
11986885	-	1257-0050	Green	50	0.8	0.24
11996885	-	1257-0100	Green	100	0.8	0.15
11906895	-	1257-0200	Blue	200	0.8	0.15
11916895	-	1257-0250	Purple	250	1.2	0.30
11926895	-	1257-0300	Purple	300	1.0	0.25
11936895	-	1257-0500	Purple	500	0.8	0.20
11946895	-	1257-1000	Purple	1,000	0.8	0.15

Pipettor tips, sterile, in VistaRak™



- Available in sterile and filtered sterile packs
- Pyrogen free
- RNase and DNase certified
- Filtered tips contain a hydrophobic filter which acts as a barrier to aqueous liquids and aerosols, protecting the pipettor and sample from trace amounts of carryover

Sterile

Catalogue No	Legacy No	Alt. No	Volume, µL	Pack type	Inner pack qty	Pack qty
11935762	-	4060-1032	25	Rack	192	980
11975762	-	4060-2132	250	Rack	192	980
11905772	-	4060-3132	1,250	Rack	192	980

Filtered, sterile

Catalogue No	Legacy No	Alt. No	Volume, µL	Pack type	Inner pack qty	Pack qty
11945762	-	4060-1332	25	Rack	192	980
11985762	-	4060-2332	250	Rack	192	980
11915772	-	4060-3332	1,250	Rack	192	980

