









Flask lids

Multi-socket/flat flange. For reaction vessels and culture vessels.



	Catalogue No	Legacy No	Alt. No	Centre socket size	Side socket size	Side socket size [5°]	Side socket size [10°]	Side socket size [15°]
1	11321875	-	DA5/75	34/35	-	-	-	-
2	10064184	-	MAF1/75	19/26	14/23	19/26	19/26	19/26

	Catalogue No	Legacy No	Alt. No	Centre socket size	Side socket size	Side socket size [5°]	Side socket size [10°]	Side socket size [15°]
3	11301875	-	DA4/100	29/32	-	-	-	-
	11311875	-	DA5/100	34/35	-	-	-	-
4	11342025	-	MAF2/2	19/26	19/26	2 x 19/26	19/26	-
5	11352025	-	MAF2/32	19/26	19/26	24/29	19/26	-
6	11362025	-	MAF2/52	19/26	19/26	19/26	19/26	34/35
7	11372025	-	MAF3/52	24/29	19/26	19/26	19/26	34/35
8	10714537	-	MAF4/41	29/32	14/23	14/23	14/23	29/32
9	11382025	-	MAF5/2	34/35	19/26	-	2 x 19/26	-



Flask support



- Ideal support for all types of round bottom flasks and vessels
- Can accept flask sizes from 50mL to 3L
- 1 size fits all - no need to buy specific sizes to match individual flasks
- Made from grey chemically resistant polymer
- Stable and safe on the bench

Catalogue No	Legacy No	Alt. No	
12787173	-	48618	

