




Flasks, round bottom, multiple neck



Round bottom, centre and one angled neck

Catalogue No	Legacy No	Alt. No	Capacity, mL	Centre socket size	Side socket size	Height, mm
12438986	-	FR25/1S/1A	25	14/23	14/23	77
12403631	-	FR50/1S/1A	50	14/23	14/23	89
12478966	-	FR100/3S/1A	100	24/29	14/23	110
12030643	-	FR100/3S/2A	100	24/29	19/26	110
12418996	-	FR250/3S/1A	250	24/29	14/23	130
12070653	-	FR250/3S/2A	250	24/29	19/26	130
12020673	-	FR500/3S/1A	500	24/29	14/23	150
12050673	-	FR500/3S/2A	500	24/29	19/26	150
11301935	-	FR1L/3S/1A	1,000	24/29	14/23	180
12070643	-	FR1L/3S/2A	1,000	24/29	19/26	180
11371945	-	FR2L/5S/2A	2,000	34/35	19/26	230

Round bottom, three parallel necks

Catalogue No	Legacy No	Alt. No	Capacity, mL	Centre socket size	Side socket size	Height, mm	Side-centre sockets distance [min.], mm
12060653	-	FR250/3S/22P	250	24/29	19/26	130	35
12040673	-	FR500/3S/22P	500	24/29	19/26	150	40
12449016	-	FR500/4S/22P	500	29/32	19/26	155	40
12060643	-	FR1L/3S/22P	1,000	24/29	19/26	180	46
11341935	-	FR1L/5S/22P	1,000	34/35	19/26	195	50
11351935	-	FR1L/5S/33P	1,000	34/35	24/29	195	50
12458996	-	FR2L/3US/22P	2,000	24/29	19/26	210	62
11391945	-	FR2L/5S/22P	2,000	34/35	19/26	230	56
11301955	-	FR2L/5S/33P	2,000	34/35	24/29	230	56
11331965	-	FR3L/5S/22P	3,000	34/35	19/26	250	65
11341965	-	FR3L/5S/33P	3,000	34/35	24/29	250	65
11391965	-	FR5L/5US/33P	5,000	34/35	24/29	250	75

Round bottom, centre and two angled side necks

Catalogue No	Legacy No	Alt. No	Capacity, mL	Centre socket size	Side socket size	Height, mm
12591738	-	FR25/1S/11A	25	14/23	14/23	77
12483621	-	FR25/2S/11A	25	19/26	14/23	80
12493621	-	FR50/1S/11A	50	14/23	14/23	89
12413631	-	FR50/2S/11A	50	19/26	14/23	92
12424501	-	FR50/2S/22A	50	19/26	19/26	92
11301945	-	FR100/2S/11A	100	19/26	14/23	110
12020643	-	FR100/3S/11A	100	24/29	14/23	110
12488966	-	FR100/3S/22A	100	24/29	19/26	110
11371955	-	FR250/2S/22A	250	19/26	19/26	130
12408996	-	FR250/3S/11A	250	24/29	14/23	130
12050653	-	FR250/3S/22A	250	24/29	19/26	130
12419016	-	FR500/3S/11A	500	24/29	14/23	150
12030673	-	FR500/3S/22A	500	24/29	19/26	150
10228054	-	FR1L/3S/22A	1,000	24/29	19/26	180
11331935	-	FR1L/5S/22A	1,000	34/35	19/26	195
12010663	-	FR2L/3US/22A	2,000	24/29	19/26	205
11381945	-	FR2L/5S/22A	2,000	34/35	19/26	230
11321965	-	FR3L/5S/22A	3,000	34/35	19/26	250

Round bottom, centre and septum side neck

Catalogue No	Legacy No	Alt. No	Capacity, mL	Centre socket size	Side socket size	Height, mm
12463621	-	FR10/1S/SEP	10	14/23	Septum	60
12473621	-	FR25/1S/SEP	25	14/23	Septum	77