


Gloves, latex, disposable, powder free, Satin Plus, Class I Medical



Length 230mm, thickness 0.13mm (palm) 0.16mm (textured fingertip). Natural colour.

| Catalogue No | Legacy No | Size | Pack qty |
|--------------|-----------|-------------|----------|
| 12425190 | - | Extra small | 100 |
| 11729843 | - | Small | 100 |
| 11739843 | - | Medium | 100 |
| 11769853 | - | Large | 100 |
| 11740663 | - | Extra large | 90 |



Gloves, neoprene, NeoTouch™, disposable, powder free, green, PPE Category III



NeoTouch™ gloves have a polyurethane inner coating that ensures easy donning, textured fingertips provide a secure grip in both dry and wet conditions. The beaded cuff ensures a secure hold, and their green colour makes them easy to identify. Highly resistant to acids, bases and alcohols. They provide good chemical splash resistance.

240mm standard length

| Catalogue No | Legacy No | Size | Pack qty |
|--------------|-----------|------|----------|
| 12475060 | - | S | 100 |
| 12485060 | - | M | 100 |
| 12495060 | - | L | 100 |
| 12405070 | - | XL | 100 |

290mm long length

| Catalogue No | Legacy No | Size | Pack qty |
|--------------|-----------|-------------|----------|
| 12415070 | - | Small | 100 |
| 12425070 | - | Medium | 100 |
| 12435070 | - | Large | 100 |
| 12445070 | - | Extra large | 100 |



Gloves, neoprene, disposable, sterile, Suretech™ Pharma, PPE Category III



- Powder free neoprene glove
- Designed for users sensitive to natural rubber latex proteins
- Packed individually in sealed pair packs
- Gamma irradiation sterilised
- Unique Biogel® coating ensures the gloves are comfortable to wear and easy to don
- Packaged in a powder free controlled environment to minimise particulate contamination
- Ideal for use in pharmaceutical manufacture
- CE Complex design

| Catalogue No | Legacy No | Size | Pack type | Pack qty |
|--------------|-----------|------|-----------|----------|
| 12910735 | - | 6 | Pair | 50 |
| 12920735 | - | 6.5 | Pair | 50 |
| 12930735 | - | 7 | Pair | 50 |
| 12940735 | - | 7.5 | Pair | 50 |
| 12950735 | - | 8 | Pair | 50 |
| 12960735 | - | 8.5 | Pair | 50 |
| 12970735 | - | 9 | Pair | 40 |

