



Disposal bins, glass, safety box, Fisherbrand®



Made from sturdy corrugated cardboard and designed for separating glass waste from other types of waste. Supplied flat packed to save storage space but can be assembled very quickly. Each box has a thick polyethylene bag, which acts as a liner for holding the glass waste and any liquid spills. The unique integral lid provides added safety. When the box is full the lid is removed and a rear flap is flipped up to seal the unit ready for disposal.

Catalogue No	Legacy No	Description	Width, mm	Depth, mm	Height, mm	Pack qty
11507692	-	Large floor model	300	300	680	6
11308503	-	Medium floor model	300	300	580	6
11318503	-	Small bench model	200	200	260	6



Disposal bin, glass, Danibin



- Orange with lid and printed label
- Easily distinguished from other waste containers
- Ideal for collection and transportation of glass and crockery
- Can be disposed of at a landfill site with other domestic refuse

Catalogue No	Legacy No	Alt. No	Capacity, L
12322449	-	DD820	22



Disposal bins, glass



Disposable, highly visible receptacle that helps reduce injuries by providing a safe laboratory work area. Made from sturdy corrugated cardboard, the receptacle has a 0.05mm thick polyethylene bag to contain broken glass fragments and spilled liquids. A pair of tabs keep the liner from slipping out of place. The snug-fitting lid has a flap through which waste material is easily inserted. When box is full, the internal safety cap is pulled into place and the entire unit is ready for disposal. The broken glass disposal box is printed with eye-catching words and graphics in four languages. Packaged flat to save space and shipping costs, it can be assembled in seconds.

Catalogue No	Legacy No	Alt. No	Description	Width, mm	Depth, mm	Height, mm	Pack qty
11302204	-	F24653-0001	Floor type	300	300	680	6
11312204	-	F24653-0002	Bench type	200	200	260	6